

Page

M1117 10-2-12 Turbidity working stds. into 20mls pipet from #M1115  
 TA (200 NTU)

6el Std=

4.00 NTU

48.9 NTU

525 NTU

2.0 mL into 18mL D.I. H<sub>2</sub>O = 20 NTU10mL into 10 mL D.I. H<sub>2</sub>O = 100 NTU20 mL into 0mL DI H<sub>2</sub>O = 200 NTU

(Exp Aug 2014).

M1118 10/02/12 LSF

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)
 1.0ml of 1000mg/L Hg (#M10563) into 100ml=10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

M1119 10/02/12 LSF

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)
 1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L into 100ml = 100ug/L Hg

M1120 10/02/12 LSF

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M 1119 ) into 100ml = 3.0ug/L Hg

M1121 10/02/12 LSF

Hg CCV  
(0.2% HNO<sub>3</sub>, HCl)
 1.0ml of 1000mg/L Hg (#M10563) into 100ml=10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

M1122 10/02/12 LSF

Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl)

0.2ml of 100ug/L Hg (#M1119 ) into 100ml = 0.2ug/L Hg

M1123 10/02/12 LSF

Hg Aqua Regia Reagent:

"Carefully" in hood, mix 3 parts HCl (AB: 553) with 1part HNO<sub>3</sub> (AB: 554)

M1124 10/04/12 LSF

Stannous Chloride Reagent:

Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (#M 1057) in 70mL HCl (AB: 553). Bring up to volume with D.I. Water.

M1125 10/04/12 LSF Potassium persulfate: 50g potassium persulfate (#M10573) dissolve into 1,000mL D.I. H<sub>2</sub>O.  
 Exp: 07/04/13

M1126 10/09/12 LSF TCEP Extraction Fluid #1: Fill a 20L Carboy w/19 L D.I. H<sub>2</sub>O, Add 114mL Glacial Acetic Acid (AR#S34) and 128.6mL 16N NaOH (#M10910). Dilute to 20L with D.I. H<sub>2</sub>O and mix pH 4.9± 0.05

M1127 10/10/12 LSF

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)
 1.0ml of 1000mg/L Hg (#M10563) into 100ml=10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

M1128 10/10/12 LSF

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)
 1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L into 100ml = 100ug/L Hg

#M11104 09/24/12 LSF TCEP Extraction Fluid #: Into a 20L carboy (inner)  
 Glacial Acetic Acid (AB:53) 128.6mL 10N NaOH (#M10563)  
 to 70L with DI H<sub>2</sub>O  
 pH = 4.93 ± 0.05

Exp: 09/24/13

#M11105 09/25/12 LSF NaCl Hydroxylamine Sulfate Reagent: Dissolve 60g NaCl (inner)  
 and 60g hydroxylamine sulfate (#M102107) in 500mL

#M11106 09/26/12 LSF Potassium Permanganate Soln: 50g Potassium permanganate (inner)  
 dissolve into 1,000mL of DI H<sub>2</sub>O

#M11107 09/26/12 LSF Hg Working Stds:  
 (0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10563) into 100ml=10mg/L Hg
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg
0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg
1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg
2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg
4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg
5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg
10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M11108 09/26/12 LSF Hg Alt Source Working Std:  
 (0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M11109 09/26/12 LSF Hg ICV / LCSW  
 (0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M11108) into 100ml = 3.0ug/L Hg
---

#M11110 09/26/12 LSF Hg CCV  
 (0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10563) into 100ml=10mg/L Hg
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg
3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M11111 09/26/12 LSF Hg MRL:  
 (0.2% HNO<sub>3</sub>, HCl)

0.2ml of 100ug/L Hg (#M11108) into 100ml = 0.2ug/L Hg
---

#M11112 09/26/12 LSF Hg Aqua Regia Reagent:

"Carefully" in hood, mix 3 parts HCl (AB: 553) with 1 part HNO<sub>3</sub> (AB: 553)

#M11113 09/26/12 LSF Stannous Chloride Reagent: Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (#M102107) in 70 mL HCl (AB: 553). Bring up to volume with D.I. Water.

#M11114 09/26/12 LSF Stannous Chloride (Macro) ~2g NO. 0000012403 Exp 11/14

#M11115 9-28-12 Stab1 Cal Standard, 200 ATU Lot A2216 ex Aug 2014

#M11116 09/28/12 LSF TCEP Extraction Fluid #: Into a 20L carboy (inner)

Acetic Acid (AB:53) 128.6mL

(#M10563) Dilute to 70L

mix. pH = 4.93 ± 0.05

#M11117

#M11118

#b (ex 10/13)	#M11094 09/19/12 (3P)	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg
17 AL (ex 11/13), .000 ug/l std	#M11095 09/19/12 (5P)	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M1092) into 100ml = 0.2ug/L Hg
97 AL (ex 2/13)			
VI10897 AL (ex 11/13)			
	#M11096 09/19/12 (5P)	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 553) with 1 part HNO <sub>3</sub> (AB: 554)
	#M11097 09/19/12 (5P)	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (#M10563) in 70 mL HCl (AB: 553). Bring up to volume with D.I. Water.
	#M11098 09-21-12 GFAA LOD/LOQ soil working std:	Into a 200ml vol. flask, partially filled with AMA (1% HNO <sub>3</sub> ) EXP. 11-25-12	D.I. H <sub>2</sub> O, pipet the following and bring up to volume with D.I. H <sub>2</sub> O. 1.0 ml Sb 1000mg/L (#M10230) = 1,000 ug/L Sb 3.2 ml As 1000mg/L (#M10281) = 10,000 ug/L As 3.2 ml Se 1000mg/L (#M10277) = 10,000 ug/L Se 0.03 ml Pb 10,000mg/L (#M10539) = 4,000 ug/L Pb 1.0 ml TL 1000mg/L (#M10566) = 5,000 ug/L TL
Ex. # M110868 d do mark w/ a DI	#M11099 09-21-12 GFAA LOD/LOQ soil Spiking soin:	Into a 1 L vol. flask, pipet 10mL of #M11098 and AMA (1% HNO <sub>3</sub> ) EXP. 11-25-12	bring up to volume with D.I. H <sub>2</sub> O.
Glacial Acetic Acid G#M10910, Dilute to H = 4.93 ± 0.05	#M11100 09-21-12 GFAA LOD/LOQ soil Ag Working std:	Into a 200ML vol. flask, partially filled with AMA (1% HNO <sub>3</sub> ) EXP. 12-3-12	D.I. H <sub>2</sub> O, pipet the following and bring up to volume with D.I. H <sub>2</sub> O. 0.12mL Ag 1000mg/L (#M10310) = 400 ug/L Ag
piped 1.0ml of H (#M10910) = 4.93 ± 0.05	#M11101 09-21-12 GFAA LOD/LOQ soil Ag Spiking soin:	Into a 1L vol. flask, pipet 10mL of #M11100 and AMA (1% HNO <sub>3</sub> ) EXP. 12-3-12	bring up to volume with D.I. H <sub>2</sub> O.
nl=10mg/L Hg lg lg lg lg lg lg - Hg	#M11102 09-21-12 GFAA LOD/LOQ water working std:	Into a 100ml vol. flask, partially filled with AMA (1% HNO <sub>3</sub> ) EXP. 11-25-12	D.I. H <sub>2</sub> O, pipet the following and bring up to volume with D.I. H <sub>2</sub> O. 1.5mL Sb 1000mg/L (#M10230) = 15,000 ug/L Sb 1.0mL As 1000mg/L (#M10281) = 10,000 ug/L As 2.0mL Se 1000mg/L (#M10277) = 20,000 ug/L Se 0.045mL Pb 10,000mg/L (#M10539) = 4500 ug/L Pb 0.4mL TL 1000mg/L (#M10566) = 4000 ug/L TL 0.1mL Ag 1000mg/L (#M10310) = 1000 ug/L Ag
l = 10mg/L Hg	#M11103 09-21-12 GFAA LOD/LOQ water Spiking soin:	Into a 1L vol. flask, pipet 10mL of #M11102 and bring AMA (1% HNO <sub>3</sub> )	up to volume with D.I. H <sub>2</sub> O.

#M11086 cont. -

100K	10 mls of #M11051 Cu (ex 2/14), 10 mls of #M11054 Mn (ex 2/14), 10 mls of #M11053 Cr (ex 2/14), 10 mls #M10839 Pb (ex 10/13)
	10 mls #M10284 Zn (ex 11/12) = 100,000 ug/L
100,000	10 mls of #M10899 Mg (ex 06/13), 10 mls of #M10896 Fe (ex 11/13), 10 mls of #M10898 Ca (ex 11/13), 10 mls #M10897 AL (ex 11/13)
	10 mls #M10570 Na (ex 4/13), 10 mls #M10570 Na (ex 3/13), 10 mls #M10657 K (ex 6/13) into 1000 ml vol flask = 100,000 ug/l std
500,000	50 mls of #M10899 Mg (ex 11/13), 50 mls of #M10896 Fe (ex 11/13), 50 mls of #M10898 Ca (ex 11/13), 50 mls #M10897 AL (ex 2/13)
	into 1000 ml vol flask = 500,000 ug/l std
1000k	100 mls of #M10899 Mg (ex 06/13), 100 mls of #M10896 Fe (ex 11/13), 100 mls of #M10898 Ca (ex 11/13), 100 mls #M10897 AL (ex 2/13)
	into 1000 ml vol flask = 1,000,000 ug/l std
	Bring the 1000 ml volumetric up to mark with DI H <sub>2</sub> O

#M11087

Continuing Calibration Standard 1(CCV1)

9-17-12 7A

Into a one liter volumetric flask pipet 50mls of # M11031 (ex 8/13)  
and 5 ml of #M11085 (ex 09/13).

ex.nov 8/13

Bring up to mark with DI H<sub>2</sub>O5% HCl  
50% HNO<sub>3</sub>

Continuing Calibration Standard 2(CCV2)

Into a one liter volumetric flask pipet 5 mls of # M11031 (ex 8/13)  
and 0.5 ml of #M11085 (ex 09-13).

ex.nov 8/13

Bring up to mark with DI H<sub>2</sub>O

#M11088 ICSA B onto 500 ml vol flask pipet 50mls #M10866  
 9-17-12 2A  
 2.5mls #M11031 + #M11085 Bring up to mark in  
 H<sub>2</sub>O ex 4/13.

\*M11089 09/17/12 LSP TCLP Extraction Fluid #1: Into a 20L carboy pipet 114ml Glacial Acetic Acid (AB-534) and 128.6ml 10N NaOH (#M10910). Dilute to 20L with D.I. H<sub>2</sub>O and mix pH = 4.93 ± 0.05  
 exp 09/17/13

\*M11090 9-18-221 TCLP Extraction fluid #1 Into 20L carboy pipet 114ml Glacial Acetic Acid (AB-534) + 128.6 ml 10N NaOH (#M10910)  
 ex 9/18/13 Dilute to 20L with DI H<sub>2</sub>O and mix pH = 4.93 ± 0.05

\*M11091 09/19/12 LSP

Hg Working Stds:  
 (0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

\*M11092 09/19/12 LSP

Hg Alt Source Working Std:  
 (0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

\*M11093 09/19/12 LSP

Hg ICV / LCSW  
 (0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M1092) into 100ml = 3.0ug/L Hg

#M11076 09/11/21SF	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10568) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg
#M11077 09/11/21SF	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg
#M11078 09/11/21SF	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M11077) into 100ml = 3.0ug/L Hg
#M11079 09/11/21SF	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg
#M11080 09/11/21SF	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M11077) into 100ml = 0.2ug/L Hg
#M11081 09/11/21SF	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 553) with 1 part HNO <sub>3</sub> (AB: 554)
#M11082 09/11/21SF	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (#M1055) in 70 mL HCl (AB: 553). Bring up to volume with D.I. Water.
#M11083 09/11/21SF	NaCl Hydroxylamine sulfate Reagent: Dissolve 60g NaCl (#M574) and 60g hydroxylamine sulfate (#M1020) in 500mL milli-Q H <sub>2</sub> O. Exp: 03/11/13	
#M11084 09/17/21SF	K <sub>2</sub> MnO <sub>4</sub> Sol'n: 50g potassium permanganate (#M9903) Exp: 03/17/13	Dissolve into 1,000mL of D.I. H <sub>2</sub> O.
#M11085 9-17-12	Cerium Assurance std Spec Corrig'd lot # XE739120 # 15-182CR ox 9-15-13 :	
#M11086 9-17-12		
<b>ICAP Working Standards</b>		
5% HNO <sub>3</sub> 5%HCL		
Standard name	Pipet the following respectively:	
0.25	0.25 mls of standard 100 ug/L into 1000 ml volumetric flask = 0.25 ug/L std	
0.5	0.50 mls of standard 100 ug/L into 1000ml volumetric flask= 0.50 ug/L std	
1	1.0 mls of standard 100ug/L into 1000 ml volumetric flask = 1.0 ug/L std	
5	5.0 mls of standard 100ug/L into 1000 ml volumetric flask =5.00 ug/L std	
10	10.0 mls of standard 100ug/L into 1000 ml volumetric flask =10.0 ug/L std	
20	20.0 mls of standard 100ug/L into 1000 ml volumetric flask =20.0 ug/L std	
50	50.0 mls of standard 100ug/L into 1000 ml volumetric flask =50.0 ug/L std	
100	1.0 mls of #M10837(ex 3/13) and 1.0 mls of #M11085 (ex 09/13) into 1000 ml volumetric flask= 100 ug/L std	
1000	10.0 mls of #M10837(ex 4/13) and #M11085 (ex 9/13) and 1 mL of #M10564 Si into 1000 ml volumetric flask= 1000 ug/L std	
10,000	100 mls #M10837 (ex 4/13), 10 mls #M10568 Be (ex 4/13), 10 mls #M10278 Cd (ex 11/12)	
6	#M11086 Cont... ↗	

300

#M11041 09/07/12 LSP

3010 LOQ Spiking Solution  
Base SPIKE PREPARATION 50

5%  $\text{HNO}_3$   
5%  $\text{HCl}$

Of this Base standard, pipet 10 mls into 500 ml volumetric to create working std or 1 ml into 50 ml digestion tube for a digested working standard.

Expires: 04/2011

M11072 09/03/2015

5% HNO<sub>3</sub>  
5% HCl

3010 LOQ Spiking Solution  
Base SPIKE PREPARATION 50

Into a 1000 mL Volumetric Flask, pipet the following:

Analyte	LOQ (ug/L)	Std ID #	Std Conc (mg/L)	Amount (mL) pipet into
K	500	M1000	10000	2.5
Na	600	M1000	10000	3.0

Of this Base standard, pipet 10 ml's into 500 ml volumetric to create working std or 1 ml into 50 ml digestion tube for a digested working standard.

Expires: 04/20/

#M11073 09/07/12 LSP

5% HNO<sub>3</sub>  
5% HCl

301B LOQ Spiking Solution  
Bass SPIKE PREPARATION 5%

Into a 1000 mL Volumetric Flask, pipet the following:

Analyte	LOQ ( $\mu\text{g/L}$ )	Std ID #	Std Conc ( $\text{mg/L}$ )	Amount ( $\text{mL}$ ) pipet into 1 L
B	46	ml0155	.1000	0.6
SI	200	ml3636	.1000	10

Of this Base standard, pipet 10 mls into 500 ml volumetric to create working std or 1 ml into 50 ml digestion tube for a digested working standard.

Expires: 11/25/17

M11624 09/10/15F - TCEP Extraction Fluid #: Into a 2.0L carboy pipette

Exp: 09/10/13

Acetic Acid (AR 534) and  
CH<sub>3</sub>CO<sub>2</sub>H). Dilute to 20 l with  
mix. pH = 4.93 ± 0.05

\* M11025 09/10/12 (3F) Stannous Chloride (JT Baker) - Lot No. 0000000579

Ex: 03/29/17

ous Chloride (#M1085) in  
ater.

#M11067 09/07/12 LSF

up to

5% HNO<sub>3</sub>  
5% HCl

3050 LOQ Spiking Solution  
Base SPIKE PREPARATION 50x

Into a 1000 mL Volumetric Flask, pipet the following:

Analyte	LOQ (ug/L)	Std ID #	Std Conc (mg/L)	Amount (mL) to pipet into 1 L
K	2040		10000	13.2
Na	960		10000	4.8

Of this Base standard, pipet 10 mL into 500 mL volumetric to create a working std or 1 mL into 50 mL digestion tube for a digested working standard.

Expires: 04/30/13

(L) to  
1 L

#M11068 2 9-7-12

ICV List 2. Into 1 Liter Volumetric flask →  
5% HCl pipet 1 mL of #M10709, #M11056, #M10567, #M10562, #M10569  
5% HNO<sub>3</sub> Bring up to mark with DI H<sub>2</sub>O. = 1000 ug/L List 2 by April 20, 2013.

date a  
rking

#M11069 09/07/12 LSF

5% HNO<sub>3</sub>  
5% HCl

3050 LOQ Spiking Solution  
Base SPIKE PREPARATION 50x

Into a 1000 mL Volumetric Flask, pipet the following:

Analyte	LOQ (ug/L)	Std ID #	Std Conc (mg/L)	Amount (mL) to pipet into 1 L
B	52	M10235	1000	2.6
Si	192	M10360	1000	9.6

Of this Base standard, pipet 10 mL into 500 mL volumetric to create a working std or 1 mL into 50 mL digestion tube for a digested working standard.

Expires: 11/25/12

(mL) to  
1 L

#M11070 09/07/12 LSF

5% HNO<sub>3</sub>  
5% HCl

3010 LOQ Spiking Solution  
Base SPIKE PREPARATION 50x

Into a 1000 mL Volumetric Flask, pipet the following:

Analyte	LOQ (ug/L)	Std ID #	Std Conc (mg/L)	Amount (mL) to pipet into 1 L
Ag	4	M10210	1000	0.2
Al	38	M10630	10000	0.38
As	24	M10781	1000	1.2
Ba	1.8	M10562	1000	0.09
Be	0.6	M10568	1000	0.03
Ca	100	M10464	10000	0.5
Cd	2	M10213	1000	0.1
Co	4	M10561	1000	0.2
Cr	4	M10553	10000	0.02
Cu	7	M11051	10000	0.035
Fe	100	M10594	10000	0.5
Mg	40	M10662	10000	0.2
Mn	4	M11054	10000	0.02
Mo	7	M10560	1000	0.35
Ni	6	M10574	1000	0.3
Pb	4	M10524	10000	0.02
Sb	12	M10230	1000	0.6
Se	13	M10211	1000	0.65
Tl	15	M10566	1000	0.75
V	5	M10565	1000	0.25
Zn	10	M10384	10000	0.05

Of this Base standard, pipet 10 mL into 500 mL volumetric to create a working std or 1 mL into 50 mL digestion tube for a digested working standard.

Expires: 11/25/12

create a  
working

**Stannous Chloride Reagent:** Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride in 70 mL HCl (AB: 553). Bring up to volume with H<sub>2</sub>O.

#M11064 09/05/12  $\text{NiNO}_3$  Modifier - Into a 100mL vol. flask, pipet 20mLs  $\text{NiNO}_3$  ( $M_{9853}$ ) and b  
Add volume with milli-Q  $\text{H}_2\text{O}$ .

#M11065 09/07/12 LSF 3056 106 ✓

**3050 LOQ Spiking Solution  
Base SPIKE PREPARATION 50**

Into a 1000 mL Volumetric Flask, pipet the following:

Analyte	LOQ ( $\mu\text{g/L}$ )	Std ID #	Std Conc (mg/L)	Amount in pipet
Ag	4	M10310	1000	0.0
Al	9.6	M10646	10000	0.0
As	32	M10781	1000	0.0
Ba	2	M11052	1000	0.0
Be	1.6	M10562	1000	0.0
Ca	56	M10494	10000	0.0
Cd	1.8	M10234	1000	0.0
Co	9.6	M10501	1000	0.0
Cr	5.6	M10553	10000	0.0
Cu	16	M10503	10000	0.0
Fe	72	M10641	10000	0.0
Mg	32	M10682	10000	0.0
Mn	6	M10104	10000	0.0
Mo	9.6	M10560	1000	0.0
Ni	4.8	M10515	1000	0.0
Pb	10	M10545	10000	0.0
Sb	32	M10280	1000	0.0
Se	16	M10617	1000	0.0
Tl	18	M10560a	1000	0.0
V	3.2	M10535	1000	0.0
Zn	12	M10502	10000	0.0

Of this Base standard, pipet 10 ml's into 500 ml volumetric flask to make working std or 1 ml into 50 ml digestion tube for a digestion standard.

Expires: 1/25/12

M10066 09/07/12 5F 3050 1,002

**3659 LOQ Spiking Solution  
Base SPIKE PREPARATION 50x**

Into a 1000 mL Volumetric Flask, pipet the following:

Of this 8-ml standard, pipet 10 mls into 500-ml volumetric flask to make working std or 1 ml into 50-ml digestion tube for oxidation standard.

Expires: 04/2011

rt prep lot #	#M11046 08/24/17 SF	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 553) with 1part HNO3 (AB: 554)
Hg	#M11047 08/24/17 SF	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (#M1062) in 70 mL HCl (AB: 553). Bring up to volume with D.I. Water.
L Hg	# M11048 08-26-12 TCCP Soln #1: Add to a 20L Cambay pipet 114 mL Glacial Acet 100% (AB: 534) and 128.6 mL 10N NaOH (#M10910). Dilute to 20L with DI H <sub>2</sub> O and mix pH = 4.9±0.0		
- Hg	#M11049 08/27/17 SF Potassium persulfate soln: 50g potassium persulfate (#M10573) dissolved into 1,000mL DI H <sub>2</sub> O. Exp: 02/27/13		
L Hg	#M11050 08/27/17 SF Potassium permanganate soln: 50g potassium permanganate (#M19903) dissolved into 1,000mL of DI H <sub>2</sub> O. Exp: 02/27/13		
- Hg	#M11051 08-28-12 CPI Cu 10,000mg/L Lot # 12G118 ex Feb 22, 2014.		
L Hg	#M11052 08-28-12 CPI Ba, 1000mg/L Lot # 12G245 ex Feb 22, 2014.		
part HNO3 (AB: 554)	#M11053 08-28-12 CPI Cr, 10,000 mg/L Lot # 11K059 ex Feb 22, 2014.		
ous Chloride (#M10985) in later.	#M11054 8-28-12 CPI Mn 10,000mg/L Lot # 11L109 ex Feb 22, 2014.		
	#M11055 8-30-12 Spex Centriprep Custom Assurance std Lot # 15-1534C ex Aug 30, 2013.		
Glacial Acetic 10N NaOH (#M10910) D and mix	#M11056 8-30-12 SCP Science Tungsten 1000ug/L ex. 12/13		
to 100ml=10mg/L Hg 0ug/L Hg 5ug/L Hg 10ug/L Hg 20ug/L Hg 50ug/L Hg 100ug/L Hg	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 50ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg	
100ml = 10mg/L Hg 0ug/L Hg	#M11057 09/04/17 SF	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L into 100ml = 100ug/L Hg
0 100ml = 3.0ug/L Hg	#M11058 09/04/17 SF	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M1058 ) into 100ml = 3.0ug/L Hg
to 100ml=10mg/L Hg 0ug/L Hg 5ug/L Hg	#M11059 09/04/17 SF	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg
0 100ml = 0.2ug/L Hg	#M11060 09/04/17 SF	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M11058 ) into 100ml = 0.2ug/L Hg
	#M11061 09/04/17 SF	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 553) with 1part HNO3 (AB: 554)
	#M11062 09/04/17 SF		

#M11031 # 8-17-12 Custom Assurance Std Spec  
S-1224P - exp Aug 15, 2013

#M11032 08/20/12 SF Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl) 1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M11033 08/20/12 SF Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl) 1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M11034 08/20/12 SF Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl) 3.0ml of 100ug/L Hg (#M 11033 ) into 100ml = 3.0ug/L Hg

#M11035 08/20/12 SF Hg CCV  
(0.2% HNO<sub>3</sub>, HCl) 1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M11036 08/20/12 SF Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl) 0.2ml of 100ug/L Hg (#M 11033 ) into 100ml = 0.2ug/L Hg

#M11037 08/20/12 SF Hg Aqua Regia Reagent: "Carefully" in hood, mix 3 parts HCl (AB: 553) with 1 part HNO<sub>3</sub> (AB: 553)

#M11038 08/20/12 SF Stannous Chloride Reagent: Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride in 70 mL HCl (AB: 553). Bring up to volume with Distilled Water.

#M11039 08/23/12 13F (Fisher) Hydroxylamine Sulfate 16+ # 11678 Exp: 08/23/13

#M11040 08/23/12 SF TCLP Extraction Fluid #1: Into a 20L carboy pipet 1/2ml Glacial Acetic Acid (AB:534) and 100ml H<sub>2</sub>O

Exp: 08/23/13  
Dilute to 20L with H<sub>2</sub>O  
pH = 2.93

#M11041 08/24/12 SF Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl) 1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M11042 08/24/12 SF Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl) 1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M11043 08/24/12 SF Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl) 3.0ml of 100ug/L Hg (#M 11033 ) into 100ml = 3.0ug/L Hg

#M11044 08/24/12 SF Hg CCV  
(0.2% HNO<sub>3</sub>, HCl) 1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M11045 08/24/12 SF Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl) 0.2ml of 100ug/L Hg (#M 11033 ) into 100ml = 0.2ug/L Hg

#M11016 08/06/12 AMA	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M11013) into 100ml = 0.2ug/L Hg
#M11017 08/06/12 AMA	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 553) with 1part HNO <sub>3</sub> (AB: 554)
#M11018 08/06/12 AMA	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (#M10985) in 70 mL HCl (AB: 553). Bring up to volume with D.I. Water.

#M11019 08/07/12 SF TCP Extraction Fluid #1: Into a 20L Carboy pipet 114ml Glacial Acetic Acid (AB:534) and 178.6ml 10N NaOH (#M10910). Dilute to 20L with D.I. H<sub>2</sub>O and mix. pH = 4.93 ± 0.05  
Exp: 08/07/13

#M11020 8-8-12 Custom Assurance std, Spox catprep Lot # 8-974P exp Aug 2013  
24

#M11021 08/09/12 SF Potassium Permanganate Sol'n: 50g potassium permanganate (#M9903) dissolve into 1,000ml of D.I. H<sub>2</sub>O.  
Exp: 02/09/13

#M11022 08/09/12 Pd/Mg Modifier - Into a 50mL vol-flask, pipet 15mL #M10726 and 5mL  
AMA  
#M9551. Bring up to volume with Milli-Q H<sub>2</sub>O. Exp. 10-30-12

#M11023 08/09/12 SF NaCl Hydroxylamine sulfate reagent: Dissolve 60g NaCl (#M8574), and 60g hydroxylamine sulfate (#M10210) in 500mL milli-Q H<sub>2</sub>O.  
Exp: 02/09/13

#M11024 08/13/12 SF	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg
---------------------	---	---

#M11025 08/13/12 SF	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg
---------------------	---	--

#M11026 08/13/12 SF	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M1025) into 100ml = 3.0ug/L Hg
---------------------	--	--

#M11027 08/13/12 SF	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg
---------------------	---	--

#M11028 08/13/12 SF	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M1025) into 100ml = 0.2ug/L Hg
---------------------	--	--

#M11029 08/13/12 SF	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 553) with 1part HNO <sub>3</sub> (AB: 554)
---------------------	------------------------	--

#M11030 08/13/12 SF	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (#M10985) in 70 mL HCl (AB: 553). Bring up to volume with D.I. Water.
---------------------	----------------------------	--

#M11002 07/23/12(5F)

Hg Aqua Regia Reagent: "Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO<sub>3</sub> (AB: 538)

#M11003 07/23/12(5F)

Stannous Chloride Reagent: Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M1095) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

#M11004 07/26/12(5F) - TCDP Extraction Fluid #1: Into a 20L carboy, pipet 114ml 60% Acetic Acid (AB: 534) and 178.6mL 10N NaOH (#10910). Dilute to 20L with deionized water. pH = 4.93 ± 0.05

Exp: 07/26/13

#M11005 07/30/12(5F)

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10663) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M11006 07/30/12(5F)

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M11007 07/30/12(5F)

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M11006) into 100ml = 3.0ug/L Hg

#M11008 07/30/12(5F)

Hg CCV  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M11009 07/30/12(5F)

Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl)

0.2ml of 100ug/L Hg (#M11006) into 100ml = 0.2ug/L Hg

#M11010 07/30/12(5F)

Hg Aqua Regia Reagent:

"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO<sub>3</sub> (AB: 538)

07/30/12(5F)

Stannous Chloride Reagent:

Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M1095) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

07/30/12(5F)

#M11011 07/30/12(5F)

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10663) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M11012 08/06/12

AMIA

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M11014 08/06/12

AMIA

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M11013) into 100ml = 3.0ug/L Hg

#M11015 08/06/12

AMIA

Hg CCV  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10663) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M10987 07/12/12 LSF Potassium Persulfate <sup>07/12/12 SF</sup> ~~salt, sog potas~~ <sup>07/12/12 SF</sup> Lot No: 117410 Exp: 07/12/13

#M10988 07/16/12 LSF

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M10989 07/16/12 LSF

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M10990 07/16/12 LSF

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M10989) into 100ml = 3.0ug/L Hg

#M10991 07/16/12 LSF

Hg CCV  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M10992 07/16/12 LSF

Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl)

0.2ml of 100ug/L Hg (#M10989) into 100ml = 0.2ug/L Hg

#M10993 07/16/12 LSF

Hg Aqua Regia Reagent:

"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO<sub>3</sub> (AB: 538)

#M10994 07/16/12 LSF

Stannous Chloride Reagent:

Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (#M10785) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

#M10995 07/17/12 LSF

TCLP Extraction Fluid #: 1: Into a 700L carboy pipet 114mL Glacial Acetic Acid (AB: 534) and 128.6mL 10N NaOH (#M10910). Dilute to 70L with D.T. H<sub>2</sub>O and mix, pH = 4.93 ± 0.05  
Exp: 07/17/13

#M10996 07/20/12 LSF

Potassium permanganate 3d'n: 30g potassium permanganate (#M9903)  
Exp: 01/20/13 dissolve into 1,000mL of D.T. H<sub>2</sub>O.

#M10997 07/20/12 LSF

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M10998 07/20/12 LSF

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M10999 07/20/12 LSF

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M10998) into 100ml = 3.0ug/L Hg

#M11000 07/20/12 LSF

Hg CCV  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

Hg MRI

0.2ml of 100ug/L Hg into 100ml = 0.2ug/L Hg

#M10976 07/09/12 15F Hg Aqua Regia Reagent: "Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO<sub>3</sub> (AB: 538)  
#M10977 07/09/12 15F Stannous Chloride Reagent: Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M 1082) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

#M10978 07/09/12 15F Potassium persulfate soln: 50g potassium persulfate (#M10573) dissolve into 1000mL D.I. H<sub>2</sub>O.  
Exp: 01/09/13

#M10979

7-10-12

2A

### ICAP Working Standards List 2

.5% HNO<sub>3</sub> .5% HCl

#### Standard

name	Pipet the following respectively:
1	0.01 mls of #M10837 and .001mls #M10559 = 1 ug/L
10	0.10 mls of #M10837 and .01mls #M10559 = 10 ug/L
100	1 mls of #M10334 and 0.1mls #M10559 = 100 ug/L
1000	10 mls of #M10334 and 1.0mls #M10559 = 1000 ug/L
10000	100 mls of #M10334 and 10mls #M10559 = 10000 ug/L

ex3-13

#M10980 7-10-12

2A

#### Initial Calibration Standard (ICV) List 2

Into 1 liter vol flask pipet 1 ml of #M10657, M10562, M10569, M10709, and M9930, bring up to mark with DI H<sub>2</sub>O ex 8/12

#M10981 7-10-12

2A

#### Continuing Calibration Standard (CCV)

into 1 liter vol flask pipet 10 mls of #M10837 and 1.0mls #M10559 = 1000 ug/L bring up to mark with DI H<sub>2</sub>O ex 3-13

#M10982 07-10-12 GFAA LOD/LOQ Soil Working Std: Into a 200mL vol. flask, partially filled with  
 AMA (1% HNO<sub>3</sub>) D.I. H<sub>2</sub>O, pipet the following and bring up to  
 EXP. 11/25/12 volume with D.I. H<sub>2</sub>O.

$$1.2 \text{ mL Sb } 1000 \text{ mg/L } (\#M10280) = 1,200 \mu\text{g/L Sb}$$

$$3.2 \text{ mL As } 1000 \text{ mg/L } (\#M10281) = 3,200 \mu\text{g/L As}$$

$$3.2 \text{ mL Se } 1000 \text{ mg/L } (\#M10277) = 3,200 \mu\text{g/L Se}$$

$$0.08 \text{ mL Pb } 10,000 \text{ mg/L } (\#M10839) = 800 \mu\text{g/L Pb}$$

$$1.0 \text{ mL Tl } 1000 \text{ mg/L } (\#M10566) = 1,000 \mu\text{g/L Tl}$$

$$0.12 \text{ mL Ag } 1000 \text{ mg/L } (\#M10310) = 120 \mu\text{g/L Ag}$$

#M10983 07-10-12 GFAA LOD/LOQ Soil Spiking soln: Into a 1L vol. flask, pipet 10ml of

AMA (1% HNO<sub>3</sub>) bring up to volume with D.I. H<sub>2</sub>O.

Exp. 11/25/12

#M10984 07/10/12 LSF Electrolyte - Lot NO.: 1W350A Exp: 12/16/14

#M10985 07/10/12 LSF Stannous chloride - Lot No. K475 88 Exp: 07/11/13

#M10986 07/10/12 LSF Nacl Hydroxylamine Sulfate Reagent: Dissolve 10g NaCl  
 hydroxylamine sulfate

#M10962 06/29/12 LSF Potassium permanganate soln: 56g potassium permanganate (#M9903) dissolve into 1,000ml of D.I. H<sub>2</sub>O.  
Exp: 12/29/12

#M10963 06/29/12 LSF TCLP Extraction Fluid #1: Add to a 20L carboy pipet 114mL Glacial Acetic Acid (AB:534) and 128.6mL 10N NaOH (#M10910). Dilute to 20L with D.I. H<sub>2</sub>O and mix. pH = 9.8 ± 0.05  
Exp: 06/29/13

of 60g  
50mL

#M10964 07/02/12 LSF

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M10965 07/02/12 LSF

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M10966 07/02/12 LSF

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M10965) into 100ml = 3.0ug/L Hg

#M10967 07/02/12 LSF

Hg CCV  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M10968 07/02/12 LSF

Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl)

0.2ml of 100ug/L Hg (#M10965) into 100ml = 0.2ug/L Hg

#M10969 07/02/12 LSF

Hg Aqua Regia Reagent:

"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO<sub>3</sub> (AB: 538)

#M10970 07/02/12 LSF

Stannous Chloride Reagent:

Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M 10241) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

#M10971 07/09/12 LSF

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M10972 07/09/12 LSF

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M10973 07/09/12 LSF

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M10972) into 100ml = 3.0ug/L Hg

#M10974 07/09/12 LSF

Hg CCV  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M10975 07/09/12 LSF

Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl)

0.2ml of 100ug/L Hg (#M10972) into 100ml = 0.2ug/L Hg

#M10954 06/25/12/5F

Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl)

0.2ml of 100ug/L Hg (#M 10951) into 100ml = 0.2ug/L Hg

#M10955 06/25/12/5F

Hg Aqua Regia Reagent: "Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO<sub>3</sub> (AB: 536)

#M10956 06/25/12/5F

Stannous Chloride Reagent: Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride in 70 mL HCl (AB: 537). Bring up to volume with DI Water.

#M10957 06/25/12/5F NaCl Hydroxylamine Sulfate Reagent: Dissolve bag NaCl hydroxylamine sulfate (AB: 537) in

Exp: 12/25/12  
Milli - o H<sub>2</sub>O.

#M10958 6-27-12

## ICAP Working Standards

5% HNO<sub>3</sub>      5% HClStandard  
name

Pipet the following respectively:

0.25	0.25 mls of standard 100 ug/L into 1000 ml volumetric flask = 0.25 ug/L std
0.5	0.50 mls of standard 100 ug/L into 1000ml volumetric flask= 0.50 ug/L std
1	1.0 mls of standard 100ug/L into 1000 ml volumetric flask =1.00 ug/L std
5	5.0 mls of standard 100ug/L into 1000 ml volumetric flask =5.00 ug/L std
10	10.0 mls of standard 100ug/L into 1000 ml volumetric flask =10.0 ug/L std
20	20.0 mls of standard 100ug/L into 1000 ml volumetric flask =20.0 ug/L std
50	50.0 mls of standard 100ug/L into 1000 ml volumetric flask=50.0 ug/L std
100	1.0 mls of #M10837(ex 3/13) and 1.0 mls of #M10470 (ex 09/12) into 1000 ml volumetric flask= 100 ug/L std
1000	10.0 mls of #M10837(ex 4/13) and 1.0 ml of #M10310 Ag, #M10568 Be, #M10278 Cd and #M10564 Si into 1000 ml volumetric flask= 1000 ug/L std
10,000	100 mls of #M10837 (ex 4/13), 10 mls #M10568 Be (ex 4/13), 10 mls #M10278 Cd (ex 11/12)
100k	10 mls of #M10144 Cu (ex 9/12) , 10 mls of #M10141 Mn (ex 9/12) , 10 mls of #M10139 Cr (ex 9/12) , 10 mls #M10284 Zn (ex 9/12)
100,000	10 mls of #M10899 Mg (ex 06/13) , 10 mls of #M10896 Fe (ex 11/13) , 10 mls of #M10898 Ca (ex 11/13) , 10 mls #M10907 Al (ex 11/13)
500,000	10 mls #M10570 Na (ex 4/13), 10 mls #M10570 Na (ex 3/13) , 10 mls #M10667 K (ex 6/13) into 1000 ml vol flask = 500,000 ug/l std
1000k	50 mls of #M10899 Mg (ex 11/13) , 50 mls of #M10896 Fe (ex 11/13) , 50 mls of #M10898 Ca (ex 11/13) , 50 mls #M10907 Al (ex 11/13) into 1000 ml vol flask = 500,000 ug/l std
	100 mls of #M10899 Mg (ex 06/13) , 100 mls of #M10896 Fe (ex 11/13) , 100 mls of #M10898 Ca (ex 11/13) , 100 mls #M10907 Al (ex 11/13) into 1000 ml vol flask = 1,000,000 ug/l std

Bring the 1000 ml volumetric up to mark with DI H<sub>2</sub>O#M10959 ICV/CSW onto 1 liter vol flask piped 10mL #M10837  
27 6-27-12 #M10568, 0.5mL #M10567. Bring up to mark with DI H<sub>2</sub>O

#M10960 CCV1, CCV2. onto 2 1 liter Volumetrics

CCV 1 → 50 mls #M10837 and 50 mls

with DI H<sub>2</sub>O

CCV 2 → 5 mls #M10837 and 0.5 mls

mark with DI H<sub>2</sub>O. ex 9/12

#M10961 CSAB onto 500 ml vol flask piped 50 mls #M10837

5% HCl 6-27-12

#M10470 Bring up to mark w/ DI H<sub>2</sub>O. a g5% HNO<sub>3</sub>

# HCl (AB:537) 538)	# M10940 06/19/12(LSF)	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg
1, dissolve 4) in 7omL HCl with D.I. H <sub>2</sub> O.	# M10941 06/19/12(LSF)	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L into 100ml = 100ug/L Hg
	# M10942 06/19/12(LSF)	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M 10941) into 100ml = 3.0ug/L Hg
	# M10943 06/19/12(LSF)	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg
	# M10944 06/19/12(LSF)	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M 10941) into 100ml = 0.2ug/L Hg
	# M10945 06/19/12(LSF)	Hg Aqua Regia Reagent	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO <sub>3</sub> (AB: 538)
	# M10946 06/19/12(LSF)	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask; dissolve 100g Stannous Chloride (# M 10941) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.
(AB: 538)	# M10947 06/19/12(LSF) Potassium Permanganate Soln: 50g potassium permanganate (#M9903) dissolve EXP: 12/19/12		into 1,000mL of D.I. H <sub>2</sub> O.
ide (# M 10941) in	# M10948 06/19/12(LSF) TCLP Extraction fluid #1: Into a 20L Carboy, pipet 114ml Glacial Acetic Acid (AB:539) and 128.6mL 10 N NaOH (#M10910). Dilute to 70L with D.I. H <sub>2</sub> O and mix. DTH = 4.73 ± 0.05		
# M10949 06/19/12	AHA	Pd/Mg Modifier - Into a 50mL vol. Flsk, pipet 15mL *M10726 and 5mL *M10551. Bring up to volume with Milli-Q water. EXP: 10-30-12	
= 100ug/L Hg 2mL = 5.0ug/L Hg	# M10950 06/19/12(LSF)	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg
100mL = 100ug/L Hg ug/L Hg 100ug/L Hg			
P ex June 15 2013	# M10951 06/19/12(LSF)	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L into 100ml = 100ug/L Hg
	# M10952 06/19/12(LSF)	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M 10951) into 100ml = 3.0ug/L Hg
	# M10953 06/19/12(LSF)	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M10928 06/04/12 LSP

Hg Aqua Regia Reagent: "Carefully" in hood, mix 3 parts HCl (AB: 537) with one part HNO<sub>3</sub> (AB: 537)

#M10929 06/04/12 LSP

Stannous Chloride Reagent: Into a 100mL volumetric flask, dissolve 100g Stannous Chloride (#M10724) (AB: 537). Bring up to volume with D.I. Water.

#M10930 06/11/12 LSP

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M10931 06/11/12 LSP

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M10932 06/11/12 LSP

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M10531) into 100ml = 3.0ug/L Hg

#M10933 06/11/12 LSP

Hg CCV  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M10934 06/11/12 LSP

Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl)

0.2ml of 100ug/L Hg (#M10531) into 100ml = 0.2ug/L Hg

#M10935 06/11/12 LSP

Hg Aqua Regia Reagent:

"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO<sub>3</sub> (AB: 537)

#M10936 06/11/12 LSP

Stannous Chloride Reagent:

Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (#M10724) (AB: 537) with 70 mL HCl (AB: 537). Bring up to volume with D.I. Water,

#M10937 06/11/12 LSP Hg 100/100 Spiking Soln(water) - 1.0mL of 1000mg/L Hg (#M10563)

$$100\text{mL} = 10\text{mg/L Hg}$$

$$1.0\text{mL of } 10\text{mg/L Hg into } 100\text{mL} = 10\text{ug/L Hg}$$

$$5.0\text{mL of } 100\text{ug/L Hg into } 100\text{mL} = 5\text{ug/L Hg}$$

#M10938 06/11/12 LSP Hg 100/100 Spiking Soln(soil) - 1.0mL of 1000mg/L Hg (#M10563) into 100mL

$$1.0\text{mL of } 10\text{mg/L Hg into } 100\text{mL} = 10\text{ug/L Hg}$$

$$1.0\text{mL of } 100\text{ug/L Hg into } 100\text{mL} = 10\text{ug/L Hg}$$

#M10939 6-18-12 Custom Assurance std Spec Cert prep lot # 7-1884P

#M10914 05/29/12 LF	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M 10913) into 100ml = 3.0ug/L Hg
#M10915 05/29/12 LF	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg
#M10916 05/29/12 LF	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M 10913) into 100ml = 0.2ug/L Hg
#M10917 05/29/12 LF	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO <sub>3</sub> (AB: 538)
#M10918 05/29/12 LF	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (#M20241) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.
#M10919 05/29/12 LF	Potassium permanganate Soln: 50g potassium permanganate (#M9903) Exp: 11/29/12	dissolve into 1000mL of D.I. H <sub>2</sub> O.
#M10920 05/29/12 LF	NaCl Hydroxylamine sulfate Reagent: Dissolve 100g NaCl (#M9574) and 100g hydroxylamine sulfate (#M10210) in 500mL milli- $\mu$ H <sub>2</sub> O. Exp: 11/29/12	
(AB: 538)		
ide (#M10241) in	#M10921 05/30/12 LF TCLP Extraction Fluid #1: Into a 20L carboy add 114ml Glacial Acetic Acid (AB: S34) and 12.84ml 10N NaOH (#M10910). Dilute to 20L with D.I. H <sub>2</sub> O and mix. pH = 4.93 ± 0.05 Exp. 05/30/13	
M9574) and (#M10210) in	#M10922 05/30/12 LF 10N NaOH: Into a 1L vol. flask add 114ml NaOH (#M10571) and bring up to volume with D.I. H <sub>2</sub> O.	
#M10923 05/30/12 LF	Hg working stds: 1.0ml of 1000mg/L Hg (#M10563) into 100mL = 10mg/L Hg (0.2% HNO <sub>3</sub> , HCl) 1.0 ml of 10mg/L Hg into 100mL = 100ug/L Hg 0.5 ml of 100ug/L Hg into 100mL = 0.5ug/L Hg 1.0 mL of 100ug/L Hg into 100mL = 1.0ug/L Hg 2.0 mL of 100ug/L Hg into 100mL = 2.0ug/L Hg 4.0 mL of 100ug/L Hg into 100mL = 4.0ug/L Hg 5.0 mL of 100ug/L Hg into 100mL = 5.0ug/L Hg 10.0 mL of 100ug/L Hg into 100mL = 10.0ug/L Hg	
#M10924 05/30/12 LF	Hg Alt Source working std: 1.0ml of 1000 mg/L Hg (#M9063) into 100mL = 10mg/L Hg (0.2% HCl/H <sub>3</sub> PO <sub>4</sub> /HCl) 1.0ml of 10mg/L Hg into 100mL = 100ug/L Hg 1.0ml of 100ug/L Hg into 100mL = 10ug/L Hg	
mg/L Hg		
#M10925 06/04/12 LF	Hg TCV 1L(SW) (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M10924) into 100mL = 3.0ug/L Hg
mg/L Hg		
#M10926 06/04/12 LF	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100mL = 10mg/L Hg 1.0ml of 10mg/L Hg into 100mL = 100ug/L Hg 3.0ml of 100ug/L Hg into 100mL = 3.0ug/L Hg
mg/L Hg		
#M10927 06/04/12 LF	Hg MRL (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M10924) into 100mL = 0.2ug/L Hg

#M10901 05/21/12 LSF

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M10902 05/21/12 LSF

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L into 100ml = 100ug/L Hg

#M10903 05/21/12 LSF

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M10902) into 100ml = 3.0ug/L Hg

#M10904 05/21/12 LSF

Hg CCV  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M10905 05/21/12 LSF

Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl)

0.2ml of 100ug/L Hg (#M10902) into 100ml = 0.2ug/L Hg

#M10906 05/21/12 LSF

Hg Aqua Regia Reagent:

"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO<sub>3</sub> (AB: 537)

#M10907 05/21/12 LSF

Stannous Chloride Reagent:

Into a 1000 mL volumetric flask, dissolve 100g Stannous chloride in 70 mL HCl (AB: 537). Bring up to volume with DI Water.

#M10908 05/21/12 LSF Marl Hydrazylamine Sulfate Reagent: Dissolve 60g N-hydrazinyl-  
 (60g hydrazylamine sulfate)  
 Sodium Nitro - a H<sub>2</sub>O

Exp: 11/21/12

#M10909 05/22/12 LSF Potassium Persulfate Soln: 50g potassium persulfate  
 Exp: 11/22/12 dissolve into 1,000 ml D.I. H<sub>2</sub>O

#M10910 05/24/12 LSF 10N NaOH: Into a 1L Vol. Flask add 400g NaOH (AB: 1054)  
 Exp: 05/24/13 bring up to volume with D.T. H<sub>2</sub>O.

#M10911 05/24/12 LSF TCLP Extraction Fluid #1: Into a 20L carboy Pipet 1L  
 Acetic Acid (AB: 534) and 178L water  
 (#M10910). Dilute to 20L total  
 and mix. pH = 4.93 ± 0.05

#M10912 05/29/12 LSF

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M10913 05/29/12 LSF

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L into 100ml = 100ug/L Hg

Manganese 10 info 1,000ml	#M10880 05-09-12 TLP Extraction Fluid #2 - Into a 20L carboy, pipet 114ml Glacial Acetic Acid (AB:534). Dilute to 20L with D.I. H <sub>2</sub> O and mix. pH = 2.93 ± 0.05. Exp. 05-09-13
#M10881 05/10/12 LSF	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)
#M10882 05/10/12 LSF	1.0ml of 1000mg/L Hg (#M40563) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg
#M10883 05/10/12 LSF	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)
#M10884 05/10/12 LSF	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg
#M10885 05/10/12 LSF	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)
#M10886 05/10/12 LSF	3.0ml of 100ug/L Hg (#M 1088) into 100ml = 3.0ug/L Hg
#M10887 05/10/12 LSF	Hg CCV (0.2% HNO <sub>3</sub> , HCl)
#M10888 05/10/12 LSF	1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg
#M10889 05/10/12 LSF	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)
#M10890 05/10/12 LSF	0.2ml of 100ug/L Hg (#M 1088) into 100ml = 0.2ug/L Hg
#M10891 05/10/12 LSF	Hg Aqua Regia Reagent: "Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO <sub>3</sub> (AB: 538)
#M10892 05/10/12 LSF	Stannous Chloride Reagent: Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M 20241) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

Glacial Acetic Acid (AB:537) and 10N NaOH (#M10855). Mix. Exp: 05/11/13

Glacial Acetic Acid (AB:537) into a 20L carboy, pipet 114ml Glacial Acetic Acid (AB:534) and 128.6ml 10N NaOH (#M10855). Dilute to 20L with D.I. H<sub>2</sub>O and mix. (pH = 4.93 ± 0.05)

#M10894 05/11/12 LSF

TCLP Extractant Fluid #1: Into a 20L carboy, pipet 114ml Glacial Acetic Acid (AB:534) and 128.6ml 10N NaOH (#M10855). Dilute to 20L with D.I. H<sub>2</sub>O and mix. (pH = 4.93 ± 0.05)

#M10895 Custom Assurance Sf

Spec. lot # 7-1104P ex. May 15 2013

5-14-12

#M10896 5-14-12 Fe 10,000mg/L CPI lot #12C111 Nov 2013

82

21

#M10897 5-14-12 Al 10,000mg/L CPI lot #11L100 ex Nov 2013

72

#M10898 5-14-12 Ca 10,000mg/L CPI lot #12D128 ex Nov 2013

72

#M10899 5-14-12 Mg 10,000mg/L CPI lot #12C187 ex Nov 2013

21

#M10900 05/18/12 LSF Potassium Permanganate soln: 10g potassium permanganate (#M9903) dissolve into 1,000mL of D.I. H<sub>2</sub>O.

Exp: 11/18/12

#M10873 05/03/12 SF Potassium permanganate solution: 50g potassium permanganate (#M10880  
(#M9903) dissolve into 1,000ml AN  
of D.T. H<sub>2</sub>O.

Exp 11/03/12 SF

#M10874 05/04/12 SF

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M10881

#M10875 05/04/12 SF

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M10888.0

#M10876 05/04/12 SF

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M10875) into 100ml = 3.0ug/L Hg

#M10889.3

#M10877 05/04/12 SF

Hg CCV  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M10890.6

#M10878 05/04/12 SF

Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl)

0.2ml of 100ug/L Hg (#M10875) into 100ml = 0.2ug/L Hg

#M10891.0

#M10879 05/04/12 SF

Hg Aqua Regia Reagent:

"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO<sub>3</sub> (AB: 538)

#M10892.6

#M10880 05/04/12 SF

Stannous Chloride Reagent:

Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (#M20241) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

#M10893

#M10881 05/04/12 SF TCLP Extraction Fluid #1: Into a 20L carboy, pipet 114ml Glacial Acetic Acid (AB534), and 128.6ml 10N NaOH (#M10855).  
Exp. 05/04/13 (pH = 4.93 ± 0.05) Dilute to 20L with D.T. H<sub>2</sub>O and mix.

#M10894

#M10882 05/04/12 SF TCLP Extraction Fluid #1: Into a 20L carboy, pipet 114ml Glacial Acetic Acid (AB534), and 128.6ml 10N NaOH (#M10855). Dilute to 20L with D.T. H<sub>2</sub>O and mix. (pH = 4.93 ± 0.05)

#M10885

#M10883 05/07/12 Stannous Chloride Crystal, Macron Chemicals Lot No. K40582  
Gee

#M10886

#M10884 05/08/12 Pd/Mg Modifier: Into a 50 mL vol. flask, pipet 15 mL #M10726 and 5 mL #M10551. Bring up to volume with MilliQ water. Exp. 10-30-12

#M10887

#M10885 05/08/12 Turbidity Standard 2.0 NTU, GFS Chemicals Lot No. C148889  
Gee, Expires: 05/13

#M10888

#M10861 04/27/12 LSP	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg
#M10862 04/27/12 LSP	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg
#M10863 04/27/12 LSP	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M10862) into 100ml = 3.0ug/L Hg
#M10864 04/27/12 LSP	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg
#M10865 04/27/12 LSP	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M10862) into 100ml = 0.2ug/L Hg
#M10866 04/27/12 LSP	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO <sub>3</sub> (AB: 538)
#M10867 04/27/12 LSP	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (#M6241) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

#M10868 04/30/12 GCE Interferents' A Spec Certi Prep Lot No. 7-324P Exp. 04-30-13

#M10869 04/30/12 GCE Telp #1 Extraction Fluid : Into a 20L carboy, pipet 114mL Glacial Acetic Acid (AB:534) and 128.6mL 10N NaOH (#M10855). Dilute to 20L with D.I. H<sub>2</sub>O and mix. pH = 4.93 ± 0.05

# M10870

5-2-12

### ICAP Working Standards Boron & Silicon

Standard name	Into 1 liter Pipet the following respectively:
10	0.1 mls of #M10837 and .01mls #M10860 = 10 ug/L
50	0.5 mls of #M10837 and .05mls #M10860 = 50 ug/L
100	1 mls of #M10837 and 0.1mls #M10860 = 100 ug/L
1000	10 mls of #M10837 and 1.0mls #M10860 = 1000 ug/L
10000	100 mls of #M10837 and 10mls #M10860 = 10000 ug/L

ex 3-13

#10871

5-2-12

### Initial Calibration Standard (ICV) B, Si

Into 1 liter vol flask pipet 1 ml of #M10275 & M10564  
bring up to mark with DI H<sub>2</sub>O ex 11/12

#10872

5-2-12

### Continuing Calibration Standard (CCV)

Into 1 liter vol flask pipet 10 mls of #M10837 and 1.0mls #M10564 = 1000 ug/L  
bring up to mark with DI H<sub>2</sub>O ex 12-12

402

4/19/12

#M10847 04/18/12/SP

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M108

#M10848 04/19/12/SP

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M108

#M10849 04/19/12/SP

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M10848) into 100ml = 3.0ug/L Hg

#M108

#M10850 04/19/12/SP

Hg CCV  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M10850

#M10851 04/19/12/SP

Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl)

0.2ml of 100ug/L Hg (#M10848) into 100ml = 0.2ug/L Hg

#M10850

#M10852 04/19/12/SP

Hg Aqua Regia Reagent:

"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO<sub>3</sub> (AB: 538)

#M10850

#M10853 04/19/12/SP

Stannous Chloride Reagent:

Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (#M10241) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

#M10850

#M10854 04/19/12/SP

Potassium Permanganate Soln: 50g potassium permanganate (#M9903)  
 Exp: 10/19/12  
 dissolve into 1000ml of D.I. H<sub>2</sub>O.

#M108

#M10855 04/24/12/SP

10N NaOH: Into a 1L vol. flask add 200g NaOH (#M10588) and bring up to volume with D.I. H<sub>2</sub>O.  
 Exp: 04/24/13

#M108

#M10856 04/24/12/SP

TCLP Extraction Fluid #1: Into 20L carboy, pipet 114ml Glacial Acetic Acid (AB:534) and 128.6ml 10N NaOH (#M10708).  
 Exp: 04/24/13 (pH=4.93±0.05) Dilute to 20L with D.I. H<sub>2</sub>O and mix.

#M10856

5.

#M10857 04/25/12

Cat

SAP Ext. Fluid Acid Mixture: Into a 200mL vol. flask, add 12g H<sub>2</sub>SO<sub>4</sub> (AB:528) and 8g Cone. HNO<sub>3</sub> (AB:551). Bring up to volume with D.I. H<sub>2</sub>O.  
 Exp: 04/25/13

#M10857

#M10858 04/26/12

Cat

SAP Extraction Fluid #1: Into a 20L carboy, add D.I. H<sub>2</sub>O to mark and adjust pH to 4.20±0.05 with #M10857.

#M10858

#M10859 04/26/12

Cat

NH<sub>2</sub> Hydroxylamine Sulfate Reagent: Dissolve 10g NaOH (#M9574) and 10g Hydroxylamine Sulfate (#M10101) in 500mL Milli-a H<sub>2</sub>O.

#M10859

#M10860 4-27-2012

Cat

Silicon 1000 mg/L FISHER lot # 116650 exp 12/13

to a lume 13	#M10833 04/05/12 TCLP Extraction Fluid #1: into a 20L carboy pipet 114ml Glacial Acetic Acid (AB: 534) and 128.6ml 10N NaOH (GM10708). Dilute to 20L with D.I. H <sub>2</sub> O and mix. pH = 4.93 ± 0.05
a Volume 13.	#M10834 04/09/12 GFAA Base Std (25 µg/L) Pipet 25mL of #M10818 into a 1L vol. flask and (GAE) (5% HCl, 2% HNO <sub>3</sub> ) bring up to volume with MilliQ H <sub>2</sub> O. Exp. off-d-B
a Volume 13.	#M10835 04/09/12 GFAA CCV Std (10 µg/L) Pipet 10mL of #M10818 into a 1L vol. flask and (GAE) (5% HCl, 2% HNO <sub>3</sub> ) bring up to volume with MilliQ H <sub>2</sub> O. Exp. off-d-B,
iss/olue	#M10836 04/09/12 GFAA Alt. Source Std (10 µg/L) Pipet 50mL of #M10819 into a 1L vol. flask (GAE) (5% HCl, 2% HNO <sub>3</sub> ) and bring up to volume with MilliQ H <sub>2</sub> O. Exp. off-d-B
-	#M10837 4-11-2012 Custom Assurance std 10µg/L per Certifprep lot # 11-189 CR ex 3/13.
-	#M10838 4-11-2012 Custom Assurance std x spike 3 per Certifprep lot # 11-188CR ex 3/13
	# M10839 4-11-12 CPI Pb 10,000 mg/L Lot# 115007 ex 10/13.
#M10840 04/12/12 SF	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)  1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100µg/L Hg 0.5ml of 100µg/L Hg into 100ml = 0.5µg/L Hg 1.0ml of 100µg/L Hg into 100ml = 1.0µg/L Hg 2.0ml of 100µg/L Hg into 100ml = 2.0µg/L Hg 4.0ml of 100µg/L Hg into 100ml = 4.0µg/L Hg 5.0ml of 100µg/L Hg into 100ml = 5.0µg/L Hg 10.0ml of 100µg/L Hg into 100ml = 10.0µg/L Hg
#M10841 04/12/12 SF	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)  1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100µg/L Hg
#M10842 04/12/12 SF	Hg/ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)  3.0ml of 100µg/L Hg (#M10841) into 100ml = 3.0µg/L Hg
) in hydroxylamine H <sub>2</sub> O	#M10843 04/12/12 SF  Hg CCV (0.2% HNO <sub>3</sub> , HCl)  1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100µg/L Hg 3.0ml of 100µg/L Hg into 100ml = 3.0µg/L Hg
	#M10844 04/12/12 SF  Hg MRL: (0.2% HNO <sub>3</sub> , HCl)  0.2ml of 100µg/L Hg (#M10841) into 100ml = 0.2µg/L Hg
#M10845 04/12/12 SF	Hg Aqua Regia Reagent:  "Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO <sub>3</sub> (AB: 538)
#M10846 04/12/12 SF 7821 3	Stannous Chloride Reagent: Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M10241) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

#M10820 04/02/12 GFAA Alt Source Std (10µg/L)  
 Cat (1% HNO<sub>3</sub>) Pipet 50ml of #M10817  
 into 1L vol flask and bring up to volume  
 with MilliQ H<sub>2</sub>O.

#M10821 04/02/12 GFAA CCV Std (10µg/L)  
 Cat (1% HNO<sub>3</sub>) Pipet 10ml of #M10818 into  
 1L vol flask and bring up to volume  
 with MilliQ H<sub>2</sub>O.

#M10822 04/02/12 GFAA Pore Std (25µg/L)  
 Cat (1% HNO<sub>3</sub>) Pipet 25ml of #M10818 into  
 1L vol flask and bring up to volume  
 with MilliQ H<sub>2</sub>O.

#M10823 04/03/12 Potassium Persulfate: 50g potassium persulfate (#M9914) into 1000mL D.I. H<sub>2</sub>O.  
 Exp: 10/03/18

#M10824 04/03/12 LSF  
 Hg Working Stds:  
 (0.2% HNO<sub>3</sub>, HCl) 1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M10825 04/03/12 LSF  
 Hg Alt Source Working Std:  
 (0.2% HNO<sub>3</sub>, HCl) 1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M10826 04/03/12 LSF  
 Hg ICV / LCSW  
 (0.2% HNO<sub>3</sub>, HCl) 3.0ml of 100ug/L Hg (#M10825) into 100ml = 3.0ug/L Hg

#M10827 04/03/12 LSF  
 Hg CCV  
 (0.2% HNO<sub>3</sub>, HCl) 1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M10828 04/03/12 LSF  
 Hg MRL:  
 (0.2% HNO<sub>3</sub>, HCl) 0.2ml of 100ug/L Hg (#M10825) into 100ml = 0.2ug/L Hg

#M10829 04/03/12 LSF  
 Hg Aqua Regia Reagent: "Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO<sub>3</sub> (AB: 538)

#M10830 04/03/12 LSF  
 Stannous Chloride Reagent: Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (AB: 536) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

#M10831 04/03/12 LSF  
 NaCl Hydroxyamine Sulfate Reagent: Dissolve 100g NaCl (#M9574) and 100g Hydroxyamine Sulfate (#M10720) in 500mL

#M10832 4-4-2012 Yttrium 10,000 Environmental Express Ltd  
 ex Ser

the following  
H<sub>2</sub>O

#M10809 03/22/12(3F)

Hg Alt Source working std: 1.0ml of 100ug/L Hg (#M10803) into 100ml = 10 ug/L Hg  
(0.1% HNO<sub>3</sub>, HCl) 1.0ml of 10mg/L into 100ml = 100ug/L Hg

1)

03/22/12(5F)

#M10810 03/22/12(5F)

Hg CV/LCSW 3.0ml of 100ug/L Hg (#M10809) into 100mL = 3.0ug/L Hg  
(0.1% HNO<sub>3</sub>, HCl)

1)

#M10811 03/22/12(5F)

Hg CCR 1.0 mL of 100ug/L Hg (#M10809) into 100mL = 10ug/L Hg  
(0.1% HNO<sub>3</sub>, HCl) 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

1)

into a  
volume#M10812 03/22/12(5F) Hg MRL: 0.2mL of 100ug/L (#M10809) into 100mL = 0.2ug/L Hg  
(0.1% HNO<sub>3</sub>, HCl)

21.0 Acid

28).

#M10813 03/22/12(5F) Hg Aqua Regia Reagent: Carefully in hood, mix 3 parts HCl (A3.537) with 1 part HNO<sub>3</sub>(A3.547)

1 with

CL

4

#M10815 03/23/12(5F) 3050/062 V

10mL HCl

10mL HNO<sub>3</sub>

metric

21H<sub>2</sub>O

to 500mL

COQ

3050-1 Spiking Solution  
Base SPIKE PREPARATION 50x

Into a 1000 mL Volumetric Flask, pipet the following:

Analyte	MDL (ug/L)	Std ID #	Std Conc (mg/L)	Amount (mL) to pipet into 1 L
Li	14.4		1000	0.72
Sn	20		1000	1
Sr	3.2		1000	0.16
Tl	8.6		1000	0.48
W	24		1000	1.2

If this Base standard, pipet 10 mL into 500 mL volumetric to create a working std at 1 mg/L into 50 mL digestion tube for a digested working standard.

EXPIRES 08/01/2013

in the D166D

2nd

#M10816 03/29/12 Inorganic Ventures CTI-SPK-1 Lot no. F2-MEB417037  
Gee Exp. 04-01-2013

9903)

#M10817 03/29/12 Inorganic Ventures CTI-GFCAL-1 Lot no. F2-MEB417036  
Gee Exp. 04-01-2013

9903)

L=10mg/L Hg

#M10818 04/02/12 GFAA Working Base Std. (1000ug/L): Pipet 10mL of #M10816 into a 1 L vol. flask & bring up to volume with MilliQ water. Exp. 04-01-13  
Gee (1% HNO<sub>3</sub>)

9

#M10819 04/02/12 GFAA Alt. Source Working Std (200ug/L): Pipet 20mL of #M10817 into a 1L vol flask and bring up to volume with MilliQ H<sub>2</sub>O. Exp. 04-01-13  
Gee (1% HNO<sub>3</sub>)

19

#M10801 03/10/12 GFAA Stability Check Working Standard: Into a 250 mL vol. flask  
 Gue 1% HNO<sub>3</sub>, Exp. 09-09-12 and bring up to volume with D.I. H<sub>2</sub>O  
 0.6mL of 1,000 mg/L As (AB.534) 03/12/11  
 0.04mL of 10,000 mg/L Pb (AB.534) 03/12/11  
 1.7mL of 1,000 mg/L Se (AB.534) 03/12/11  
 0.9mL of 1,000 mg/L Sb (AB.534) 03/12/11  
 0.4mL of 1,000 mg/L Tl (AB.534) 03/12/11  
 0.07mL of 1,000 mg/L Ag (AB.534)

#M10802 03/12/12 GFAA Stability Check Standard: Pipet 5mL of #M10801  
 Gue 1% HNO<sub>3</sub>, Exp. 09-09-12 1 L vol. flask and bring up to volume with D.I. H<sub>2</sub>O

#M10803 03/13/12 TCLP Extraction Fluid #: Into a 70L carboy pipet 114mL Glacial Acetic Acid (AB.534) and 128.6mL 10N NaOH (AB.534) Dilute to 70L with D.I. H<sub>2</sub>O and mix. (pH=4.93±0.05) Exp. 03/13/13

#M10804 03/13/12 1N HCl for TCLP: Into a 1L vol. flask, partially filled with MilliQ-H<sub>2</sub>O, add 83 mL of concentrated HCl (AB.543) and bring up to volume with MilliQ water.

#M10805 3-20-2012 LOD, LOQ check 5% H<sub>2</sub>O onto 1 liter volumetric flask pipet 1mL M10762 bring up to mark with MilliQ water X50 dilution - use 1mL into 50 mL 10mL volumetric flask = 200ug/L.

#M10806 3-20-2012 LOD, LOQ check (3050) onto 1 liter volumetric flask pipet 1mL M10762. Bring up to mark with MilliQ water X50 dilution - use 1mL into 50 mL 10mL volumetric flask = 240 mg/L or 6mg/kg.

#M10807 03/22/2012 Potassium permanganate std: 50g potassium permanganate (AB.534) dissolve into 1000 mL of D.I. H<sub>2</sub>O

#M10808 03/22/2012 Hg working stds:  
 (0.7% HNO<sub>3</sub>, HCl)  
 1mL of 1000 mg/L Hg (#M1083) into 100mL volumetric flask  
 1mL of 10mg/L Hg into 100mL = 100 mg/L Hg  
 0.5mL of 10mg/L Hg into 100mL = 50 mg/L Hg  
 1mL of 10mg/L Hg into 100mL = 10 mg/L Hg  
 2.0mL of 10mg/L Hg into 100mL = 20 mg/L Hg  
 4.0mL of 10mg/L Hg into 100mL = 40 mg/L Hg  
 5.0mL of 10mg/L Hg into 100mL = 50 mg/L Hg  
 10.0mL of 10mg/L Hg into 100mL = 100 mg/L Hg

#M10790 03/06/12LF	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M 10789) into 100ml = 3.0ug/L Hg
#M10791 03/06/12LF	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg
#M10792 03/06/12LF	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M 10789) into 100ml = 0.2ug/L Hg
#M10793 03/06/12LF	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO <sub>3</sub> (AB: 538)
	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M 10572) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.
#M10794 03/09/12LF	Hg Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg
#M10795 03/09/12LF	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L into 100ml = 100ug/L Hg
#M10796 03/09/12LF	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M 10795) into 100ml = 3.0ug/L Hg
#M10797 03/09/12LF	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg
#M10798 03/09/12LF	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M 10795) into 100ml = 0.2ug/L Hg
#M10799 03/09/12LF	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO <sub>3</sub> (AB: 538)
#M10800 03/09/12LF	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M 10572) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

21  
3-12-12

#M10781 02/28/12(5F)

Hg Alt Source Working Std: 1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
(0.2% HNO<sub>3</sub>, HCl)  
1.0ml of 10mg/L into 100ml = 100ug/L Hg

#M10790

#M10782 02/28/12(5P)

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)  
3.0ml of 100ug/L Hg (#M10777) into 100ml = 3.0ug/L Hg

#M10791 02/28/12(3P)

#M10783 02/28/12(3P)

Hg CCV  
(0.2% HNO<sub>3</sub>, HCl)  
1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M10792

#M10780 02/28/12(5F)

Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl)  
0.2ml of 100ug/L Hg (#M10777) into 100ml = 0.2ug/L Hg

#M10793 0

#M10781 02/28/12(5P)

Hg Aqua Regia Reagent: "Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO<sub>3</sub> (AB: 538)

#M10782 02/28/12(5P)

Stannous Chloride Reagent: Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (#M1052) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

#M10794

#M10783 03/01/12

G6

TCP Extraction Fluid #1 : Into a 20L carboy, pipet 114mL Glacial Acetic Acid (AB: 534) and 128.6 mL 10N NaOH (#M10708). Dilute to 20L with D.I. H<sub>2</sub>O and mix. Exp. 03/02/13

#M10784

3-2-2012

Sulfur Working Std ICap onto 4 100 ml vial flasks, pipet

0.4CCV respectively from #M10773 10,000 mg/L Sulfur Std.

#M10795

also CCV< 1000ug = .01mls of #M10773 = 1000 ug/L

#M10796 0

10,000ug = .1mls of #M10773 10,000ug/L

#M10797 0

100,000ug=1.0mls of #M10773 100,000 ug/L

1000K std= 10mls of #M10773. 1000K ug/L

Bringing up to mark with DI H<sub>2</sub>O. exp. 10/13

#M10798 0

#M10785

3-2-2012

Sulfur ICN - onto 100 ml Vial flask pipet .01mls

#M10798 0

#M10767 10,000 mg/L Sulfur = 1000ug/L S.

#M10799 0

#M10786 03/05/12

G6

Pd/Mg Modifier : Into a 50 ml vol. flask, pipet 15mL #M10726 and 5mL

#M10800 0

#M10551, Brng up to volume with MilliQ water. Exp. 10-30-12

#M10787 03/05/12

G6

TCP Extraction Fluid #2 : Into a 20L carboy, pipet 114mL Glacial Acetic

#M10801

Acid (AB: 534). Dilute to 20L with D.I. H<sub>2</sub>O.

and mix pH = 2.88 ± 0.05. Exp. 03-05-13

#M10788 03/06/12(3C)

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl).

1.0ml of 1000mg/L Hg (#M9063) into 100ml=10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

3/2/13

#M10789 03/06/12(5F)

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl).

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M10764 02/17/12 LSF NaCl Hydroxylamine sulfate Reagent: Dissolve 60g NaCl (#M49574) and  
60g Hydroxylamine sulfate (#M10520) in 500 mL milli-<sup>o</sup>H<sub>2</sub>O.

Exp: 08/17/12 LSF

#M10765 02/21/12 LSF	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg
#M10766 02/21/12 LSF	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg
#M10767 02/21/12 LSF	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M10766) into 100ml = 3.0ug/L Hg
#M10768 02/21/12 LSF	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg
#M10769 02/21/12 LSF	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M 10766) into 100ml = 0.2ug/L Hg
#M10770 02/21/12 LSF	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO <sub>3</sub> (AB: 538)
#M10771 02/21/12 LSF	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M 10572) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.
#M10772 02/21/12 LSF Potassium permanganate: 50g potassium permanganate (#M9948) dissolve into 1000mL D.I. H <sub>2</sub> O.	Exp: 08/21/12	
#M10773 2-21-2012	Sulfur: 10,000ug/ml	SCP Science Lot # S120116017. ex 10/13
#M10774 02/26/12	TCLP Extraction Fluid #1: Into a 20L Carboy, pipet 114ml Glacial Acetic Acid (AB534) and 128.6ml 10N NaOH (#M10708). Dilute to 20L with D.I. H <sub>2</sub> O and mix. pH = 4.93 ± 0.05	
#M10775 02/27/12	TCLP Extraction Fluid #1: Into a 20 L Carboy, pipet 114 mL Glacial Acetic Acid (AB534) and 128.6 mL 10N NaOH (#M10708). Dilute to 20 L with D.I. H <sub>2</sub> O and mix. pH = 4.93 ± 0.05	
#M10776 02/28/12	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg
1/2/12 13S	t/C Acid Dilute to 05'	

#M10748 02/07/12/5F	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L into 100ml = 100ug/L Hg	#M107
#M10749 02/07/12/5F	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M 10748) into 100ml = 3.0ug/L Hg	#M10765C
#M10750 02/07/12/5F	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg	#M10765C
#M10751 02/07/12/5F	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M 10748) into 100ml = 0.2ug/L Hg	
#M10752 02/07/12/5F	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO <sub>3</sub> (AB: 538)	#M10766
#M10753 02/07/12/5F	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M6592) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.	
#M10754 02/14/12/5F	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg	#M10767
#M10755 02/14/12/5F	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L into 100ml = 100ug/L Hg	#M1077C
#M10756 02/14/12/5F	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M 10755) into 100ml = 3.0ug/L Hg	#M107
#M10757 02/14/12/5F	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg	#M107
#M10758 02/14/12/5F	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M 10755) into 100ml = 0.2ug/L Hg	#M107
#M10759 02/14/12/5F	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO <sub>3</sub> (AB: 538)	#M107
#M10760 02/14/12/5F	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M6592) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.	#M107 13F
#M10761 2-15-12	Phosphorus 10,000 mg/L	Spxx Centriprep (cat # AC12-161P)	#M1077 G16
#M10762 2-15-12	Sulfur 10,000 mg/L	Spxx Centriprep (cat # AC12-163S)	#M1077
#M10763 02/16/12/5F	TCLp Extractant (cat # 1): Into a 20mL Carboy Pipet 114mL Glacial Acetic Acid (AB: 534) and 128.6mL 10 N NaOH (#M 10538) Dilute to 20L with DI H <sub>2</sub> O and mix. pH = 4.93 ± 0.05		#M1077
Exp: 02/16/13			

#M10734 01-23-12 GSE	Hg CCV: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg
#M10735 01-23-12 GSE	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M10732) into 100ml = 0.2ug/L Hg
#M10736 01-23-12 GSE	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO <sub>3</sub> (AB: 538)
#M10737 01-23-12 GSE	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M10572) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.
3326 (a) in	#M10738 01-23-12 GSE	Potassium Permanganate: 50g Potassium Permanganate (# M9903) dissolved into 1,000 mL P.E. H <sub>2</sub> O. Exp. 07-23-12
SCR	#M10739 01/23/12SF TCLP Extraction Fluid #1:	To a 20L carboy pipet 114ml Glacial Acetic Acid (AB539) and 12.86mL 10N NaOH (#M10708). Dilute to 20L with D.I. H <sub>2</sub> O and mix. pH = 4.93 ± 0.05
#M10740 01/23/12LSF	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg
10726	#M10741 01/23/12SF	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)
Flask, 12.	#M10742 01/23/12SF	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)
1. Flst. 04-01-12	#M10743 01/23/12SF	Hg CCV (0.2% HNO <sub>3</sub> , HCl)
Flask, 01-12	#M10744 01/23/12SF	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)
#M10745 01/23/12SF	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO <sub>3</sub> (AB: 538)
#M10746 01/23/12SF	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M10572) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.
#10747 01/23/12SF	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

2/2

#M10720 01/17/12 SF

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M 10719) into 100ml = 3.0ug/L Hg

#M10721 01/17/12 SF

Hg GCV  
(0.2% HNO<sub>3</sub>, HCl)1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M10722 01/17/12 SF

Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl)

0.2ml of 100ug/L Hg (#M 10719) into 100ml = 0.2ug/L Hg

#M10723 01/17/12 SF

Hg Aqua Regia Reagent:

"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO<sub>3</sub> (AB: 536)

#M10724 01/17/12 SF

Stannous Chloride Reagent:

Into a 1000 mL volumetric flask, dissolve 100g-Stannous Chloride (# M10326 (a)) in 70 ml HCl (AB: 537). Bring up to volume with D.I. Water.

# M10725 1-17-12 2a Spent Cert prep Cadmium Assurance Std. lot # 10-170CR  
Gve 1/15/2013#M10726 01-17-12 Matrix Modifier 1. Pd in 8% (wt/wt) HNO<sub>3</sub> + Tr. HCl

Gve

Exp. Exp. Lot No. 1101221 Exp. 01-12-13

#M10727 01-17-12 Pd/Mg Matrix Modifier: Into a 50ml vol. flask, pipet 15ml #M10726

Gve

and 5ml #M10551. Bring up to volume with  
MilliQ water. Exp. 10-30-12

#M10728

#M10728 Gve GFAA Base Std (25ug/L) Pipet 12.5ml of #M10184 into a 250 ml vol. flask.  
01/20/12 (17. HNO<sub>3</sub>) Bring up to volume with MilliQ H<sub>2</sub>O. Exp. 04-01-12

#M10729

#M10729 01-20-12 GFAA A4 Source Std (25ug/L) Pipet 31.25 mL of #M10185 into a 250 mL vol. flask.  
Gve (17. HNO<sub>3</sub>) Bring up to volume with MilliQ H<sub>2</sub>O. Exp. 04-01-12

#M10730

#M10730 01-20-12 GFAA C4V Std (25ug/L) Pipet 50ml of #M10728 into a 100 mL vol. flask.  
Gve (17. HNO<sub>3</sub>) Bring up to volume with MilliQ H<sub>2</sub>O. Exp. 04-01-12

#M10731

#M10731 01-23-12 Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)1.0ml of 1000mg/L Hg (#M10563) into 100ml=10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M10732

#M10732 01-23-12 Gve

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M10733

#M10733 01-23-12 Gve

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M10732) into 100ml = 3.0ug/L Hg

#M 10707 01/06/12 TCLP extraction Fluid #1: Into a 20L carboy, add 114ml glacial Acetic Acid (AB: 537) and 178.6mL 10N NaOH (#M10538). Dilute to 20L with D.I. H<sub>2</sub>O and mix. pH = 4.93 ± 0.05

#M10708 01/06/12 10N NaOH: Into a 1L vol. flask add 400g NaOH (#M10537) and bring up to volume with D.I. H<sub>2</sub>O

#M10709 01/06/12 Titanium 1000 µg/mL CPI Lot No. 11E121 Exp. 04-20-13  
GCE Rec'd 10-25-11

#M10710 01/10/12 (SF)

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10668) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100µg/L Hg  
0.5ml of 100µg/L Hg into 100ml = 0.5µg/L Hg  
1.0ml of 100µg/L Hg into 100ml = 1.0µg/L Hg  
2.0ml of 100µg/L Hg into 100ml = 2.0µg/L Hg  
4.0ml of 100µg/L Hg into 100ml = 4.0µg/L Hg  
5.0ml of 100µg/L Hg into 100ml = 5.0µg/L Hg  
10.0ml of 100µg/L Hg into 100ml = 10.0µg/L Hg

labeled

#M10711 01/10/12 (SF)

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100µg/L Hg

100

#M10712 01/10/12 (SF)

Just (mL) to  
get into 1 L  
2.5  
2.5

Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100µg/L Hg (#M 6711 ) into 100ml = 3.0µg/L Hg

#M10713 01/10/12 (SF)

Hg QCV  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9984) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100µg/L Hg  
3.0ml of 100µg/L Hg into 100ml = 3.0µg/L Hg

do in to

#M10714 01/10/12 (SF)

Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl)

0.2ml of 100µg/L Hg (#M 6711 ) into 100ml = 0.2µg/L Hg

part of the

#M10715 01/10/12 (SF)

Hg Aqua Regia Reagent:

"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO<sub>3</sub> (AB: 536)

#M10716 01/10/12 (SF)

Stannous Chloride Reagent:

Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M10328 (a)) in 70 mL H<sub>2</sub>O (AB: 537). Bring up to volume with DI Water.

10.0mL

#M10717 01/16/12

Custom Assurance Standard 100 µg/L As, As<sub>3</sub>, Ba, Ca, Co, Cr, Cu, Fe, Li, Hg, Mn, Mo, Ni, Pb, Sb, Se, Sn, Sr, Ti, Tl, V, Zn. Spec Certi Prep Lot No. X-CTW1-5-500 5-157YP  
01-16-12 GCE

Exp 01-15-13

flask

#M10718 01/13/12 (SF)

Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M10668) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100µg/L Hg  
0.5ml of 100µg/L Hg into 100ml = 0.5µg/L Hg  
1.0ml of 100µg/L Hg into 100ml = 1.0µg/L Hg  
2.0ml of 100µg/L Hg into 100ml = 2.0µg/L Hg  
4.0ml of 100µg/L Hg into 100ml = 4.0µg/L Hg  
5.0ml of 100µg/L Hg into 100ml = 5.0µg/L Hg  
10.0ml of 100µg/L Hg into 100ml = 10.0µg/L Hg

cm

#M10719 01/13/12 (SF)

Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100µg/L Hg

#M10703 1-4-12 MDL Spiking Soln' all .50X

into 1 Liter val flask Pipet the following:

Analyte	(ug/L)	Std ID #	Std Conc (mg/L)	Amount (mL) to pipet into 1 L
Ag	5	M10310	1000	0.25
Al	10	M10870	10000	0.05
As	20	M10281	1000	1
Ba	2	M10143	1000	0.1
Be	0.5	M10568	1000	0.025
Ca	50	M10669	10000	0.25
Cd	0.5	M10278	1000	0.025
Co	2.5	M10581	1000	0.125
Cr	5	M10139	10000	0.025
Cu	10	M10144	10000	0.05
Fe	20	M10671	10000	0.1
Mg	25	M10688	10000	0.125
Mn	5	M10141	10000	0.025
Mo	2.5	M10560	1000	0.125
Ni	5	M10585	1000	0.25
Pb	10	M10142	10000	0.05
Sb	20	M10280	1000	1
Se	20	M10277	1000	1
Tl	20	M10568	1000	1
V	5	M10285	1000	0.25
Zn	5	M10284	10000	0.025

ex- 9-9-12

Bring up to mark with DI H<sub>2</sub>O. use 1ml into 50 mls or 100 mls  
ex. 9-9-12#M10704 1-4-12 mol Spiking Soln' Na, K. into One Liter Val  
flask pipet the following

Analyte	(ug/L)	Std ID #	Std Conc (mg/L)	Amount (mL) to pipet into 1 L
K	500	M10667	10000	1
Na	500	M10570	10000	1

Bring up to mark with DI H<sub>2</sub>O. use 1ml into bomb or 10.0 mls  
ex. 4-2-2013

#M10705 1-4-12 MDL Spiking Soln B, S: into 1 Liter flask pipet the following

Analyte	(ug/L)	Std ID #	Std Conc (mg/L)	Amount (mL) to pipet into 1 L
B	20	M10225	1000	1
S	100	M10584	1000	5

Bring up to mark with DI H<sub>2</sub>O. use 1ml into bomb or 10.0 mls  
ex. NOV 25, 2013#10706 1-4-12 m DC Spiking Soln. lists into 1 Liter val flask  
pipet the following

Analyte	(ug/L)	Std ID #	Std Conc (mg/L)	Amount (mL) to pipet into 1 L
Li	5	M10587	1000	0.25
Sn	5	M10582	1000	0.25
Sr	5	M10583	1000	0.25
Tl	5	M10580	1000	0.25
W	10	M10589	1000	0.5

Bring up to mark with DI H<sub>2</sub>O. use 1ml into 50 mls or 100 mls  
ex. 4-20-2013

#M10679 12/19/11 Gee	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M10676) into 100ml = 3.0ug/L Hg
#M10678 12/19/11 Gee	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg
#M10679 12/19/11 Gee	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M10676) into 100ml = 0.2ug/L Hg
#M10680 12/19/11 Gee	Hg Aqua Regia Reagent	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO <sub>3</sub> (AB: 538)
#M10681 12/19/11 Gee	Stannous Chloride Reagent	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M10326 (a)) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.
#M10682 12/20/11 Gee	GFAA LOD/LOD Soil Working Std (1% HNO <sub>3</sub> ) Expires: 09/09/12	Into a 200 ml vol. flask, partially filled with H <sub>2</sub> O, pipet the following and bring up to volume with D.I. H <sub>2</sub> O. 1.2 mL Sb 1000 mg/L (#M10280) = 6,000 ug/L Sb 3.2 mL As 1000 mg/L (#M10281) = 16,000 ug/L As 3.2 mL Se 1000 ug/L (#M10277) = 16,000 ug/L Se 0.08 mL Pb 10,000 mg/L (#M10141) = 4,000 ug/L Pb 1.0 mL Ti 1,000 mg/L (#M10566) = 5,000 ug/L Ti 0.12 mL Ag 1000 mg/L (#M10310) = 600 ug/L Ag
#M10683 12/20/11 Gee	GFAA LOD/LOD Soil Spiking Solution (1% HNO <sub>3</sub> ) Expires: 09/09/12	Into a 1 L vol. flask, pipet 10mL of #M10682 and bring up to volume with D.T. H <sub>2</sub> O.
#M10684 12/23/11 25L	Potassium Permanganate: 50g potassium permanganate (#M9203) dissolved into 1,000ml D.T. H <sub>2</sub> O. Exp:06/02/12	
#M10685 12/23/11SF	Hg Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg
#M10686 12/23/11SF	Hg Alt Source Working Std:	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg
#M10687 12/23/11SF	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M10676) into 100ml = 3.0ug/L Hg
#M10688 12/23/11SF	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml=10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M10672 12/16/11(3)

3010 LOG

39-10 LOQ Spiking Solution  
Base SPIKE PREPARATION

5% Et<sub>2</sub>NH

Into a 1000 mL Volumetric Flask, pipet the following:

GF/RB Base standard, pipet 10 mls into 500 ml volumetric to create a working std or 1 ml into 50 ml digestion tube for a digested working standard.

Expires: 07/09/11

#M10673 12/16/11 13:

30102009v

**3010 LAQ Spiking Solution  
Base SPIKE PREPARATION 50%**

5% HNO<sub>3</sub>  
5% HCl

Into a 1000 mL Volumetric Flask, pipet the following:

Analyte	LOQ ( $\mu\text{g/L}$ )	Std ID #	Std Conc ( $\text{mg/L}$ )	Amount (mL) to pipet into 1 L
K	500	M1038H	10000	2.5
Na	600	AA1057D	10000	3

Of this Base standard, pipet 10 ml's into 500 ml volumetric to create a working std or 1 ml into 50 ml digestion tube for a digested working standard.

卷之三

#M10674 12/16/11 L38

30102062

3010 LOQ Spiking Solution  
Base SPIKE PREPARATION 50

5% HNO<sub>3</sub>

Into a 1000 mL Volumetric Flask, pipet the following:

Analyte	LOQ ( $\mu\text{g/L}$ )	Std ID #	Std Conc ( $\text{mg/L}$ )	Amount (mL) to pipet into 1 L
B	16	M10235	1000	0.8
Si	200	M10236	1000	10

Of this Base standard, pipet 10 mls into 800 ml volumetric to create a working std or 1 ml into 50 ml digestion tube for a digested working standard.

Exhibit 04

110675 12/19/11  
P745

Hg Working Stds  
(0.2% HNO<sub>3</sub>, HC)

1.0ml of 1000mg/L Hg (#M10563) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

4106710 12/19/11  
66

Hg Alt Source Working Std  
(0.2% HNO<sub>3</sub>; HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M10666 12/15/11 SE 30102009 v

3010 LOQ Spiking Solution  
Base SPIKE PREPARATION 50x

#M106

Into a 1000 mL Volumetric Flask, pipet the following:

5% HNO<sub>3</sub>  
5% HCl

Analyte	LOQ (ug/L)	Std ID #	Std Conc (mg/L)	Amount (mL) to pipet into 1 L
Ag	4	M10310	1000	0.2
Al	36	M10310	10000	0.18
As	24	M10281	1000	1.2
Ba	1.8	M10143	1000	0.09
Be	0.6	M10563	1000	0.03
Ca	100	M10275	10000	0.5
Cd	2	M10273	1000	0.1
Zn	4	M10561	1000	0.2
Cr	4	M10139	10000	0.02
Cu	7	M10194	10000	0.035
Fe	100	M10414	10000	0.5
Mg	40	M10418	10000	0.2
Mn	4	M10441	10000	0.02
Mo	7	M10560	1000	0.35
Ni	6	M10565	1000	0.3
Pb	4	M10142	10000	0.02
Sb	12	M10280	1000	0.6
Se	13	M10277	1000	0.65
Tl	15	M10266	1000	0.75
V	5	M10285	1000	0.25
Zn	10	M10284	10000	0.05

Of this Base standard, pipet 10 mls into 500 ml volumetric to create a working std or 1 ml into 50 ml digestion tube for a digested working standard.

Expires: 09/09/12

#M10667 12/15/11 Potassium 10,000 ± 30 µg/mL CPI Lot No. 11J155 Expires: 06-09-13  
GSE

#M106

#M10668 12/15/11 Magnesium 10,000 ± 30 µg/mL CPI Lot No. 11K178 Expires: 06-09-13  
GSE#M10669 12/15/11 Calcium 10,000 ± 30 µg/mL CPI Lot No. 11K228 Expires: 06-09-13  
GSE#M10670 12/15/11 Aluminum 10,000 ± 30 µg/mL CPI Lot No. 11J015 Expires: 06-09-13  
GSE#M10671 12/15/11 Iron 10,000 ± 30 µg/mL CPI Lot No. 11L009 Expires: 06-09-13  
GSE

#M10663 12/15/11 LS

305010Q✓

5% HNO<sub>3</sub>  
5% HCl

**3050 LOQ Spiking Solution  
Base SPIKE PREPARATION 50**

Into a 1000 mL Volumetric Flask, pipet the following:

Of this Base standard, pipet 10 mls into 500 ml volumetric to create Working std or 1 ml into 50 ml digestion tube for a digested working standard.

Expires: 08/11

#10664 12/15/11 LTP

30501002

5% HNO<sub>3</sub>

**3050 LOQ Spiking Solution  
Base SPIKE PREPARATION 50%**

Into a 1000 mL Volumetric Flask, pipet the following:

Analyte	LOQ ( $\mu\text{g/L}$ )	Std ID #	Std Conc ( $\text{mg/L}$ )	Amount ( $\text{mL}$ ) to pipet into 1 L
K	2640	M10574	10000	13.2
Na	980	M10570	10000	4.8

Of this Base standard, pipet 10 mls into 500 ml volumetric to create a working std or 1 ml into 50 ml digestion tube for a digested working standard.

Expires: 01/23/

~~22 M10605 12/15/11(3P)~~

3050 LOAN

5% HNO<sub>3</sub>  
5% HCl

3839 LOQ Spiking Solution  
Base SPIKE PREPARATION 5

Into a 1000 mL Volumetric Flask, pipet the following:

Analyte	LOQ ( $\mu\text{g/L}$ )	Std ID #	Std Conc. ( $\text{mL/L}$ )	Amount ( $\text{mL}$ ) to pipet into 1 L
B	52	AM0125	1000	2.6
SI	182	AM0125	1000	9.6

Of this Base standard, pipet 10 ml's into 500 ml volumetric to create working std or 1 ml into 50 ml digestion tube for a digested working standard.

Expires: 04/16

#M10626 11-21-11 Gee	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg	11/21/11 1639
#M10627 11-21-11 Gee	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M10624) into 100ml = 0.2ug/L Hg	11/21/11 1640 11
#M10628 11-21-11 Gee	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO <sub>3</sub> (AB: 538)	
#M10629 11-21-11 Gee	Stannous Chloride Reagent	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M10326 (a)) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.	11/21/11 1641 12
#M10630 11-28-11 Gee (1% HNO <sub>3</sub> ) Exp. 09-09-12	GFAA LOD/LOQ Water Spike Working Standard : Into a 100 mL vol. flask, partially filled with D.I. H <sub>2</sub> O, pipet the following and bring up to volume with D.I. H <sub>2</sub> O.	1.5mL Sb (1,000 mg/L) #M10280 = 15,000 µg/L Sb 1.0mL As (1,000 mg/L) #M10281 = 10,000 µg/L As 2.0mL Se (1,000 mg/L) #M10277 = 20,000 µg/L Se 0.045mL Pb (1,000 mg/L) #M10142 = 4500 µg/L Pb 0.4mL Ti (1,000 mg/L) #M10566 = 4000 µg/L Ti 1.0mL Ag (1,000 mg/L) #M10310 = 10,000 µg/L Ag 0.1mL Cd (1,000 mg/L) #M10310 = 1,000 µg/L Cd	11/28/11 1702 11/28/11 1703 11/28/11 1704
#M10631 11-28-11 Gee (1% HNO <sub>3</sub> ) Exp. 09-09-12	GFAA LOD/LOQ Water Spiking Solution: Into a 1L vol. flask, pipet 10mL of #M10630 and bring up to volume with D.I. H <sub>2</sub> O.		11/28/11 1705
#M10632 11/29/11 15F	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg	11/29/11 1706
#M10633 11/29/11 15F	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg	11/29/11 1707
#M10634 11/29/11 15F	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M10623) into 100ml = 3.0ug/L Hg	
#M10635 11/29/11 15F	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg	
#M10636 11/29/11 15F	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M10623) into 100ml = 0.2ug/L Hg	
#M10637 11/29/11 15F	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO <sub>3</sub> (AB: 538)	
#M10638 11/29/11 15F	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M10326 (a)) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.	

#M10617

ICSA B NAK. onto 500 ml vol flask pipet 50 mls #M10550  
 mL 11-16-11 A 15mls #M10417, 5mls #M10570 + 5mls. #M10387.  
 2 12/12 .5% HNO<sub>3</sub> Bring up to mark w/ DI H<sub>2</sub>O. exp 10/12

#M10618

Working Standards for Sodium and Potassium; into seven, 200 ml volumetric flasks  
 11-16-11  
 .5% HCl & HNO<sub>3</sub>

For the 0.5 mg/L std - 0.1 mls of each

1.0 mg/L std - 0.2 mls of each

5.0 mg/L std - 1.0 mls of each

10.0 mg/L std - 2.0 mls of each

50.0 mg/L std - 10.0 mls of each

100 mg/L std - 20.0 mls each also used for Continuing Calibration Standard (CCV)

200 mg/L std - 40.0 mls of each

Bring up to mark with DI H<sub>2</sub>O. expires December 2012.

#M10619

Initial Calibration Standard (ICV) Into 500 ml volumetric pipette 5 mls of #M10570 and  
 #M10387 and bring up to mark with DI H<sub>2</sub>O. expires January 2013.

11-16-11

#M10620 ICSA B. Lot 2. onto 500 ml vol flask pipet 50 mls  
 11-16-11 A #M10550 15mls M10417 0.2mls of #M10616,  
 5mls HNO<sub>3</sub> + HCl #M10562 #M10567 #M10559, #M10569 Bring up  
 to Mark with DI H<sub>2</sub>O.

#M10621 Gc Pd/Mg Matrix Modifier: Into a 50ml vol. flask, pipet 15ml #M10032  
 11-17-11 and 5ml #M10551. Bring up to volume with  
 MilliQ H<sub>2</sub>O. Exp. 01-06-12

#M10622 Spec Custom Assurance std Spec Certiprep lot # 10-19 CR.  
 11-18-11 A

#M10623 Gc

11-21-11

Hg Working Stds:  
 (0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

#M10624 Gc

11-21-11

Hg Alt Source Working Std:  
 (0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

10417

DI #M10626 Gc

Hg ICV / LCSW  
 (0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M10624) into 100ml = 3.0ug/L Hg

# M10608 MRL Working std.

11-16-11 21

.5% HNO<sub>3</sub>+HCl

→ pipet from #10607 10 mls into 500ml

vol flask. Bring up to mark with DI H<sub>2</sub>O ex 6-12

M10608

ASAP

10/12

# M10609

11-16-11 21

.5% HNO<sub>3</sub> .5%HCL

## Standard

name	Pipet the following respectively: Into 1 liter vol flask
1	.01 mls of #M10334 and .001mls #M10559 = 1 ug/L
10	.10 mls of #M10334 and .01mls #M10559 = 10 ug/L
100	1 mls of #M10334 and 0.1mls #M10559 = 100 ug/L
1000	10 mls of #M10334 and 1.0mls #M10559 = 1000 ug/L
10000	100 mls of #M10334 and 10mls #M10559 = 10000 ug/L

ex 6/12

# M10610

11-16-11 21

## Initial Calibration Standard (ICV) List 2

Into 1 liter vol flask pipet 1 ml of #M10559, M10562, M10567, M10566, and M9930  
bring up to mark with DI H<sub>2</sub>O ex 6-12

# M10611

11-16-11 21

## Continuing Calibration Standard (CCV)

Into 1 liter vol flask pipet 10 mls of #M10334 and 1.0mls #M10559 = 1000 ug/L  
bring up to mark with DI H<sub>2</sub>O ex 6-12

# M10612

11-16-11

## ICAP Working Standards Boron &amp; Silicon

## Standard

name	Into 1 liter Pipet the following respectively:
10	0.1 mls of #M10334 and .01mls #M99150 = 10 ug/L
50	0.5 mls of #M10334 and .05mls #M99150 = 50 ug/L
100	1 mls of #M10334 and 0.1mls #M99150 = 100 ug/L
1000	10 mls of #M10334 and 1.0mls #M99150 = 1000 ug/L
10000	100 mls of #M10334 and 10mls #M99150 = 10000 ug/L

ex 6/12

# M10613

11-16-11 21

## Initial Calibration Standard (ICV) B, Si

Into 1 liter vol flask pipet 1 ml of #M10564 & M10575  
bring up to mark with DI H<sub>2</sub>O ex 6-12

11-16-11

# M10614

11-16-11 21

## Continuing Calibration Standard (CCV)

into 1 liter vol flask pipet 10 mls of #M10334 and 1.0mls #M99150 = 1000 ug/L  
bring up to mark with DI H<sub>2</sub>O ex 6-12

# M10615

11-16-11

ICSA I into 500 ml vol flask pipet 50 mls

#M10550 and 15 mls #M10447 BEING UP TO MARK

.5% HNO<sub>3</sub>+HCl DI H<sub>2</sub>O ex 10/12

# M10616

11-16-11

ICSA B into 500 ml vol flask pipet 50 mls #M10550 + 15 mls #M10447

.5% HNO<sub>3</sub>+HCl DI H<sub>2</sub>O ex 10/12

2.5 mls #M10334 or 2.5 mls M10447

Bring up to mark with

DI H<sub>2</sub>O ex 6/12

0.5 mls #M99150

110604

11-16-11

2A

## Continuing Calibration Standard 1(CCV1)

5% HNO<sub>3</sub>, 5% HCl

Into a one liter volumetric flask pipet 50mls of # M10334 (ex 6/12) and 5 ml of #M10310 Ag (ex 12/12), #M10568 Be 4/13) and #M10278 Cd (ex 11/12). ex.nov 6/12  
Bring up to mark with DI H<sub>2</sub>O

110605

11-16-11

2A

## Continuing Calibration Standard 2(CCV2)

5% HNO<sub>3</sub> + HCl

Into a one liter volumetric flask pipet 5 mls of # M10334 (ex 6/12) and 0.5 ml of #M10310 Ag (ex 12/12), #M10568 Be 4/13) and #M10278 Cd (ex 11/12). ex.nov 6/12  
Bring up to mark with DI H<sub>2</sub>O

#M10606

11-16-11

2A

## Initial Calibration Standard (ICV)

Into one liter volumetric flask pipet 10 mls of M10503 (ex 9/12) 2.0 mls of #m10550 (ex 10/12) and 0.5 mls #M10560 Mo (ex 4/13) and 0.5 mls of #M10275 B (ex 11/12) bring up to mark with DI H<sub>2</sub>O ex 6/12

3: 536)

( # M10326 (a)) in

#M10607

11-16-11

2A

pet 114ml

#M10538 5% HNO<sub>3</sub>

mix

HCl

= 1000 ug/L std

2/13)

MRL Std  
Base SPIKE PREPARATION 50x

Into a 1000 mL Volumetric Flask, pipet the following:

Analyte	MRL (ug/L)	Std ID #	Std Conc (mg/L)	Amount (mL) to pipet into 1 L	Expiration Date
Ag	10	M10310	1000	0.5	12/08/12
Al	400	M10416	10000	2	02/10/13
Ba	10	M10143	1000	0.5	09/09/12
Be	4	M10568	1000	0.2	04/20/13
Cd	5	M10278	1000	0.25	11/25/12
Co	10	M10561	1000	0.5	04/20/13
Cr	10	M10139	10000	0.05	09/09/12
Cu	10	M10144	10000	0.05	09/09/12
Mg	500	M10418	10000	2.5	02/10/13
Mn	10	M10141	10000	0.05	09/09/12
Mo	10	M10560	1000	0.5	04/20/13
Ni	10	M10565	1000	0.5	04/20/13
Pb	10	M10142	10000	0.05	09/09/12
Sb	20	M10280	1000	1	11/25/12
V	10	M10265	1000	0.5	11/25/12
Zn	10	M10284	10000	0.05	11/25/12
K	5000	M10387	10000	25	01/22/13
Na	1000	M10570	10000	5	04/20/13
As	10	M10281	1000	0.5	11/25/12
Ca	500	M10415	10000	2.5	02/10/13
Fe	300	M10417	10000	1.5	02/10/13
Se	10	M10277	1000	0.5	11/25/12
Tl	20	M10566	1000	1	04/20/13
Sr	10	M10569	1000	0.5	04/20/13
Ti	10	M10566	1000	0.5	04/20/13
Li	20	M10567	1000	1	04/20/13
Sn	10	M10562	1000	0.5	04/20/13
B	10	M10275	1000	0.5	11/25/12
Si	30	M10564	1000	1.5	04/20/13
W	50	M10559	1000	2.5	04/20/13

Of this Base standard, pipet 10 mls into 500 ml volumetric to create a working std or 1 ml into 50 ml digestion tube for a digested working standard

Bring up to mark w/DI H<sub>2</sub>O

#M10595 11/11/11SF	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg	604 6-11 24
#M10596 11/11/11SF	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg	05 11e-11 24
#M10597 11/11/11SF	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M10596) into 100ml = 3.0ug/L Hg	24
#M10598 11/11/11SF	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg	6006 11 24
#M10599 11/11/11SF	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M10596) into 100ml = 0.2ug/L Hg	24
#M10600 11/11/11SF	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO <sub>3</sub> (AB: 566)	6007
#M10601 11/11/11SF	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.	11 24

#M10602 11/14/11: TCLP Extraction Fluid #1 onto 20 L Carboy  
 Glacial Acetic Acid # AB. 534, + 128.6 mls  
 on 11-14-12 10N NaOH Dilute to 20L with DI H<sub>2</sub>O  
 pH = 4.93 ± 0.5.

#M10603 24 11-16-11 Do D Cal Stds

#### ICAP Working Standards

.5% HNO<sub>3</sub>   .5%HCL

Standard name.	Pipet the following respectively:	ex 6/12
0.25	0.25 mls of standard 100 ug/L into 1000 ml volumetric flask = 0.25 ug/L std	
0.5	0.50 mls of standard 100 ug/L into 1000ml volumetric flask= 0.50 ug/L std	
1	1.0 mls of standard 100ug/L into 1000 ml volumetric flask =1.0 ug/L std	
5	5.0 mls of standard 100ug/L into 1000 ml volumetric flask =5.00 ug/L std	
10	10.0 mls of standard 100ug/L into 1000 ml volumetric flask =10.0 ug/L std	
20	20.0 mls of standard 100ug/L into 1000 ml volumetric flask =20.0 ug/L std	
50	50.0 mls of standard 100ug/L into 1000 ml volumetric flask =50.0 ug/L std	
100	1.0 mls of #M10334(ex 6/12) and 1.0 mls of #M10470 (ex 09/12) into 1000 ml volumetric flask= 100 ug/L std	
1000	10.0 mls of #M10334(ex 6/12) and 1.0 ml of #M10310 Ag, #M10568 Be, #M10278 Cd and #M10564 Si into 1000 ml volumetric flask=1000 ug/L std	
10,000	100 mls #M10334 (ex 6/12),10 mls #M10568 Be (ex 4/13),10 mls #M10278 Cd (ex 11/12)	
100k	10 mls of #M10144 Cu (ex 9/12) , 10 mls of #M10141 Mn (ex 9/12) , 10 mls of #M10139 Cr (ex 9/12) , 10 mls #M10142 Pb (ex 9/12)	
100,000	10 mls #M10284 Zn (ex 11/12) = 100,000 ug/L	
500,000	10 mls of #M10418 Mg (ex 02/13) , 10 mls of #M10417 Fe (ex 02/13) ,10 mls of #M10415 Ca (ex 2/13),10 mls #M10416 Al (ex 2/13)	
1000k	10 mls #M10570 Na (ex 4/13) into 1000 ml vol flask = 100,000 ug/l std	
	50 mls of #M10418 Mg (ex 4/13) , 50 mls of #M10417 Fe (ex 02/13) ,50 mls of #M10415 Ca (ex 2/13), 50 mls #M10416 AL(ex 2/13)	
	into 1000 ml vol flask = 500,000 ug/l std	
	100 mls of #M10418 Mg. (ex 02/13) , 100 mls of #M10417 Fe (ex 2/13) ,100 mls of #M10415 Ca (ex 2/13),100 mls #M10416 (AB) (ex 2/13)	
	into 1000 ml vol flask = 1,000,000 ug/l std	
	Bring the 1000 ml volumetric up to mark with DI H <sub>2</sub> O	

# M10565 10-28-11 Nickel 1000mg/L CPI lot # HGD15 ex 4/13  
21

# M10566 10-28-11 Thallium 1000 mg/L CPI lot # 11E194 ex 4/13  
22

# M10567 10-28-11 Lithium 1000mg/L CPI lot # 11E088 ex 4/13  
23

# M10568 10-28-11 Beryllium 1000 mg/L CPI lot # 11F239 ex 4/13  
24

# M10569 10-28-11 Strontium 1000 mg/L CPI lot # 11D191 4/13  
25

# M10570 10-28-11 Na 10,000 mg/L CPI Lot # 111032 ex 4/13  
26

# M10571 10-28-11 NaOH Fisher Lot # 112665 ex 10-28-13  
27

# M10572 Stannous Chloride Fisher lot # 113509 ex 10-28-13  
28 10-28-11

# M10573 10-28-11 Potassium Persulfate lot # 11346 ex 10-28-13  
29

# M10574 11/01/113F Potassium permanganate: 50g Potassium permanganate (#M9903) dissolved into 1000mL D.I. H<sub>2</sub>O. Exp: 05/01/12

# M10575 11/01/113F Hg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)  
1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

# M10576 11/01/113F Hg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)  
1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

# M10577 11/01/113F Hg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)  
3.0ml of 100ug/L Hg (#M 10576) into 100ml = 3.0ug/L Hg

# M10578 11/01/113P Hg CCV  
(0.2% HNO<sub>3</sub>, HCl)  
1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg  
1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

# M10579 11/01/113F Hg MRL:  
(0.2% HNO<sub>3</sub>, HCl)  
0.2ml of 100ug/L Hg (#M 10576) into 100ml = 0.2ug/L Hg

# M10580 11/01/113P Hg Aqua Regia Reagent:  
"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1 part HNO<sub>3</sub> (AB: 536)

# M10591 11/01/113P Stannous Chloride Reagent:  
Into a 1000-mL volumetric flask, dissolve 100g Stannous Chloride (# M10026) into 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

#M10549 TCEP Extraction Fluid #1. Into 20 liter Carboy, pipet 114 ml  
 29 Glacial Acetic Acid (AB: 534) and 128.6 mls 10N NaOH #M10538  
 10-26-11 Dilute to 20L with DI H<sub>2</sub>O & mix. pH = 4.93 ± 0.05.

#M10550 - Inv-A. Spare Centprep (wt/wt Hg) ex Oct 3, 2012  
 29 10-26-11

wt

10N  
11x

#M10551 Magnesium W/trad Matrix Modifier Sln. 20,000 mg/L Hg  
 29 10-26-11 Spec Confirms Lot # 1-169 MG 2.

I bring

#M10552 10/26/11  
GCEHg Working Stds:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9954) into 100ml=10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg  
 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg  
 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg  
 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg  
 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg  
 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg

id

2e

#M10553 10/26/11  
GCEHg Alt Source Working Std:  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg

#M10554 10/26/11  
GCEHg ICV / LCSW  
(0.2% HNO<sub>3</sub>, HCl)

3.0ml of 100ug/L Hg (#M10553 ) into 100ml = 3.0ug/L Hg

#M10555 10/26/11  
GCEHg CCV  
(0.2% HNO<sub>3</sub>, HCl)

1.0ml of 1000mg/L Hg (#M9954) into 100ml=10mg/L Hg  
 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg  
 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg

#M10556 10/26/11  
GCEHg MRL:  
(0.2% HNO<sub>3</sub>, HCl)

0.2ml of 100ug/L Hg (#M10553 ) into 100ml = 0.2ug/L Hg

#M10557 10/26/11  
GCE

Hg Aqua Regia Reagent:

"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO<sub>3</sub> (AB: 536)#M10558 10/26/11  
GCE

Stannous Chloride Reagent:

Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M10326 (a)) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

#M10559 10-28-11 Tungsten 1000 mg/L CPI lot# 11E122. ex 4/13  
 29

#M10560 10-28-11 Molybdenum 1000 mg/L CPI lot# 10L090 ex 4/13  
 29

#M10561 10-28-11 Cobalt 1000 mg/L CPI lot# 11H211 ex 4/13  
 29

#M10562 10-28-11 Tin 1000 mg/L CPI lot# 11E120 ex 4/13  
 29

#M10563 10-28-11 Hg 1000 mg/L CPI lot# 11J054 ex 4/13  
 29

#M10564 10-28-11 Si 1000 mg/L CPI lot# 11F032 ex 4/13  
 29

326 (a)) in

	#M10493 09/23/11 L3F	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 0.5ml of 100ug/L Hg into 100ml = 0.5ug/L Hg 1.0ml of 100ug/L Hg into 100ml = 1.0ug/L Hg 2.0ml of 100ug/L Hg into 100ml = 2.0ug/L Hg 4.0ml of 100ug/L Hg into 100ml = 4.0ug/L Hg 5.0ml of 100ug/L Hg into 100ml = 5.0ug/L Hg 10.0ml of 100ug/L Hg into 100ml = 10.0ug/L Hg
	4ml		
	L		
2.7 H <sub>2</sub> O	#M10494 09/23/11 L3P	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg
142.0	#M10495 09/23/11 L3P	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100ug/L Hg (#M 10494) into 100ml = 3.0ug/L Hg
	#M10496 09/23/11 L3P	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100ug/L Hg 3.0ml of 100ug/L Hg into 100ml = 3.0ug/L Hg
	#M10497 09/23/11 L3P	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100ug/L Hg (#M 10494) into 100ml = 0.2ug/L Hg
	#M10498 09/23/11 L3P	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO <sub>3</sub> (AB: 536)
	#M10499 09/23/11 L3P	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M10326 (a)) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.

#M10500 09/23/11 Turbidity Working Standards : Take 20 mL's pipet from #M10200 (200 NTU)  
 GCE Gel Ex Stds.  
 3.88  
 50.7.  
 521  
 (Exp. Sep -12)

#M10501 09/26/11 NaCl Hydroxylamine Sulfate Reagent: Dissolve 60g NaCl (#M9574(C)) and 60g Hydroxylamine Sulfate (#M10210.) in 500mL MilliQ H<sub>2</sub>O. Exp. 03-26-12

#M10502 09/28/11 Potassium Permanganate - 50g Potassium Permanganate (#M9903)  
 ANA dissolved into 1000ml of DI H<sub>2</sub>O. Exp. 03/28/12

#M10503 9-30-11 Customer Assurance std Spex Certiprep lot # 9-84CR  
 exp Sept 30 2012.

#M10504 9-30-11 (Fisher) Hydroxylamine Sul fate lot # 111789 exp 9-30-12.

#M10505 10/03/11 TCEP Extraction Fluid #1 - Fill a 20L carboy with 19L DI. H<sub>2</sub>O. Add 114ml Glacial Acetic Acid (AB-534) and 128.4ml 10N NaOH (#M10290). Dilute to 20L with DI-H<sub>2</sub>O and mix. The pH = 4.93 ± 0.05. Exp. 10/03/12

#M10469 09/07/11 X SPIKE - 1-250 Spec Cert:Prep Lot No. 4-42YP Exp. 09-15-2012  
GIE

#M10470 09/07/11 X CTW1-4-500 Spec Cert:Prep Lot No. 4-41YP Exp. 09-15-2012  
GIE

#M10471 09/08/11 GFAA LOD/LOQ Water Spiking Sol'n: Into a 1 L vol. flask, pipet 10mL of  
GIE (1% HNO<sub>3</sub>) #M10286 and 1mL #M10472  
Expires 11-15-11

#M10472 09/08/11 GFAA Ag Working Std: Into a 100 mL vol. flask, partially filled with D.I.  
GIE (1% HNO<sub>3</sub>) H<sub>2</sub>O, pipet the following and bring up to volume  
Exp. 12-08-12 with D.I. H<sub>2</sub>O.  
1.0ml Ag (1000 mg/L) #M10310 = 10,000 µg/L Ag

#M10473 09/08/11 GFAA LOD/LOQ Soil Spiking Sol'n: Into a 1 L vol. flask, pipet 10mL of  
GIE (1% HNO<sub>3</sub>) #M10289 and 0.6mL #M10472. Bring  
Exp. 11-15-11 up to volume with D.I. H<sub>2</sub>O.

#M10474 09/12/11 AMA	Hg Working Stds: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100µg/L Hg 0.5ml of 100µg/L Hg into 100ml = 0.5µg/L Hg 1.0ml of 100µg/L Hg into 100ml = 1.0µg/L Hg 2.0ml of 100µg/L Hg into 100ml = 2.0µg/L Hg 4.0ml of 100µg/L Hg into 100ml = 4.0µg/L Hg 5.0ml of 100µg/L Hg into 100ml = 5.0µg/L Hg 10.0ml of 100µg/L Hg into 100ml = 10.0µg/L Hg
-------------------------	---	--

#M10475 09/12/11 AMA	Hg Alt Source Working Std: (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9063) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100µg/L Hg
-------------------------	---	--

#M10476 09/12/11 AMA	Hg ICV / LCSW (0.2% HNO <sub>3</sub> , HCl)	3.0ml of 100µg/L Hg (#M10475) into 100ml = 3.0µg/L Hg
-------------------------	--	---

#M10477 09/12/11 AMA	Hg CCV (0.2% HNO <sub>3</sub> , HCl)	1.0ml of 1000mg/L Hg (#M9954) into 100ml = 10mg/L Hg 1.0ml of 10mg/L Hg into 100ml = 100µg/L Hg 3.0ml of 100µg/L Hg into 100ml = 3.0µg/L Hg
-------------------------	---	---

20 #M10478 09/12/11 AMA	Hg MRL: (0.2% HNO <sub>3</sub> , HCl)	0.2ml of 100µg/L Hg (#M10475) into 100ml = 0.2µg/L Hg
-------------------------------	--	---

20ml #M10479 09/12/11 AMA	Hg Aqua Regia Reagent:	"Carefully" in hood, mix 3 parts HCl (AB: 537) with 1part HNO <sub>3</sub> (AB: 536)
---------------------------------	------------------------	--

2 #M10480 09/12/11 AMA	Stannous Chloride Reagent:	Into a 1000 mL volumetric flask, dissolve 100g Stannous Chloride (# M10326 (a)) in 70 mL HCl (AB: 537). Bring up to volume with D.I. Water.
------------------------------	----------------------------	--

2 #M10481 09/12/11 AMA	Hg LOD/LOQ spiking soln (water) - 1.0mL of 1000mg/L Hg (#M9954) into 100mL = 10mg/L Hg 1.0mL of 10mg/L Hg into 100mL = 100µg/L Hg 5.0mL of 100µg/L Hg into 100mL = 5.0µg/L Hg
------------------------------	---

#M10461

9-6-11 7A

## ICAP Working Standards

5% HNO<sub>3</sub>

5%HCL

## Standard

name Pipet the following respectively:

0.25	0.25 mls of standard 100 ug/L into 1000 ml volumetric flask = 0.25 ug/l. std
0.5	0.50 mls of standard 100 ug/L into 1000ml volumetric flask= 0.50 ug/l. std
1	1.0 mls of standard 100ug/L into 1000 ml volumetric flask =1.00 ug/l. std
.5	5.0 mls of standard 100ug/L into 1000 ml volumetric flask =5.00 ug/L std
10	10.0 mls of standard 100ug/L into 1000 ml volumetric flask =10.0 ug/L std
20	20.0 mls of standard 100ug/L into 1000 ml volumetric flask =20.0 ug/L std
50	50.0 mls of standard 100ug/L into 1000 ml volumetric flask =50.0 ug/L std
100	1.0 mls of #M10334(ex 6/12) and 1.0 mls of #M9816 (ex 8/11) into 1000 ml volumetric flask= 100 ug/L std
500	5.0 mls of #M10334(ex 6/12) and 5.0 mls of #M9916 (ex 8/11) into 1000 ml volumetric flask= 500 ug/L std <i># 9-6-11</i>
1000	10.0 mls of #M10334(ex 6/12) and 1.0 ml of #M10310 Ag, #M9930 Be, #M10278 Cd and #M9915 Si into 1000 ml volumetric flask= 1000 ug/L std
10,000	100 mls #M10334 (ex 6/12), 10 mls #M9930 Be ex (11/15/11), 10 mls #M10278 Cd ex 11/25/11
100k	10 mls of #M10144 (Cu) (ex 9/12) , 10 mls of #M10141 (Mn)(ex 9/12) , 10 mls of #M10139 (Cr)(ex 9/12), 10 mls #M10142 (Pb)(ex 9/12)
	10 mls #M10284 (Zn)(ex 11/12) = 100,000 ug/L
100,000	10 mls of #M10312 (Mg) (ex 12/12) , 10 mls of #M10417 (Fe) (ex 02/13) ,10 mls of #M10415 (Ca)(ex 2/13),10 mls #M10416 (Al)(ex 2/13)
	10 mls #M9932 Na ex 11/15/11 into 1000 ml vol flask = 100,000 ug/l std
500,000	50 mls of #M10312 (Mg) (ex 12/12) , 50 mls of #M10417 (Fe) (ex 02/13) ,50 mls of #M10415 (Ca)(ex 2/13),50 mls #M10416 (Al)(ex 2/13)
	into 1000 ml vol flask = 500,000 ug/l std
1000k	100 mls of #M10312 (Mg) (ex 12/12) , 100 mls of #M10417 (Fe) (ex 2/13) ,100 mls of #M10415 (Ca)(ex 2/13),100 mls #M10416 (Al)(ex 2/13)
	into 1000 ml vol flask = 1,000,000 ug/l std

Bring the 1000 ml volumetric up to mark with DI H<sub>2</sub>O

#M10462 NaCl Hydroxylamine Sulfato Reagent: Dissolve 10g NaCl (#M9574) and 60g  
*09/06/11 BP* Exp: 03/06/12 Hydroxylamine Sulfate (#M0210) in 500ml DI H<sub>2</sub>O

#M10463 09/06/11/5K Potassium persulfate: 50g Potassium Persulfate (#M9814) Dissolve in  
*09/06/11/5K* Exp: 03/06/12 1000 ml D.I. H<sub>2</sub>O.

#M10465 09/07/11 GFAA Base Std (25 µg/L) Pipet 5.0 mL of 1000 µg/L Sb (#M10184) into a 200mL  
*GIE* vol. flask and bring up to volume with M1110 H<sub>2</sub>O  
 18% HNO<sub>3</sub>  
 6% HCl Exp. 04-01-12

#M10466 09/07/11 GFAA Alt Source (10 µg/L) Pipet 10mL of 200 µg/L Sb (#M10185) into a 200mL  
*GIE* vol. flask and bring up to volume with M1110 H<sub>2</sub>O  
 18% HNO<sub>3</sub>  
 6% HCl Exp. 04-01-12.

#M10467 09/07/11 GFAA CCV Std. (10 µg/L) Pipet 2mL of 1000 µg/L Sb (#M10184) into a 200mL  
*GIE* vol. flask and bring up to volume with M1110 H<sub>2</sub>O. Exp. 04-01-12

#M10468 09/07/11 GFAA Sb MDL ✓ Pipet 4mL of 10 µg/L Sb (#M10187) into a 50mL  
*GIE* vol. flask and bring up to volume with M1110 H<sub>2</sub>O. Exp. 04-01-12